



MAR 14, 2024

OPEN ACCESS



DOI:

[dx.doi.org/10.17504/protocols.io.dm6gp34p1vzp/v1](https://dx.doi.org/10.17504/protocols.io.dm6gp34p1vzp/v1)

**Protocol Citation:** Ramhari Kumbhar, Hanseok Ko, Valina L. Dawson, Ted Dawson, Xiaobo Mao 2024. Stereotaxic  $\alpha$ -Syn PFF injection into mouse striatum. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.dm6gp34p1vzp/v1>

**License:** This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working  
We use this protocol and it's working

## Stereotaxic $\alpha$ -Syn PFF injection into mouse striatum

 Forked from [Stereotaxic  \$\alpha\$ -Syn PFF injection into mouse striatum](#)

Ramhari Kumbhar<sup>1,2,3</sup>, Hanseok Ko<sup>1,2,3</sup>, Valina L. Dawson<sup>4,5,3,6,7,8</sup>,  
Ted Dawson<sup>4,5,3,6,7,8,9</sup>, Xiaobo Mao<sup>4,5,3,6,7,8,9</sup>

<sup>1</sup>Neuroregeneration and Stem Cell Programs, Institute for Cell Engineering, The Johns Hopkins University School of Medicine, Baltimore, Maryland, USA;

<sup>2</sup>Department of Neurology, The Johns Hopkins University School of Medicine, Baltimore, Maryland, USA;

<sup>3</sup>Aligning Science Across Parkinson's (ASAP) Collaborative Research Network, Chevy Chase, MD 20815, USA;

<sup>4</sup>Neuroregeneration and Stem Cell Programs, Institute for Cell Engineering, Johns Hopkins University School of Medicine, Baltimore, MD 21205, USA;

<sup>5</sup>Department of Neurology, Johns Hopkins University School of Medicine, Baltimore, MD 21205, USA;

<sup>6</sup>Adrienne Helis Malvin Medical Research Foundation, New Orleans, LA 70130-2685, USA;

<sup>7</sup>Department of Physiology, Johns Hopkins University School of Medicine, Baltimore, MD 21205 USA;

<sup>8</sup>Solomon H. Snyder Department of Neuroscience, Johns Hopkins University School of Medicine, Baltimore, MD 21205, USA;

<sup>9</sup>Institute for NanoBioTechnology, Johns Hopkins University, Baltimore, MD, USA

ASAP Collaborative Research Network

Kaplitt Protocols



Eileen Ruth Torres  
Weill Cornell Medicine

### ABSTRACT

For stereotaxic injection of  $\alpha$ -Syn PFF into the mouse striatum using 2-3 month old mice. This forked protocol has slight variations in anesthesia solution and injection coordinations.

**Created:** Mar 14, 2024

**Last Modified:** Mar 14, 2024

**PROTOCOL integer ID:** 96713

**Keywords:** ASAPCRN

**Funders Acknowledgement:**  
Aligning Science for Parkinson's  
Grant ID: 020608

## MATERIALS

- Anesthetic- isoflurane OR pentobarbital (45 mg/kg) diluted 1:1 with 0.9% NaCl OR ketamine cocktail (ketamine 100mg/kg; xylazine 20mg/kg; acepromazine 3mg/kg i.p. or s.c.)
- Small electric shaver
- scalpel
- 26 gauge needle and Hamilton syringe
- minisyringe pump (Model BS-800, Braintree Scientific)
- sterile wound healing clips (7 mm stainless steel wound clip; Reflex 7 skin closure system) or sutures
- bupivacaine (0.4-0.8 ml/kg of a 0.25% solution) OR 0.5% xylocaine 1.4 ml/kg, injected sc, injected) OR 0.5% xylocaine

## BEFORE START INSTRUCTIONS

Personnel should use aseptic techniques and conduct the surgery using sterile gloves in a sterile environment. Surgical instruments should be sterilized before surgeries using an autoclave and during the procedure using glass bead sterilizer.

- 1 Anesthetize mice with a ketamine cocktail (ketamine 100mg/kg; xylazine 10mg/kg; i.p. or s.c.).
- 2 Shave the animals' skulls.
- 3 Create a 1.5 cm incision in the scalp after cleaning the scalp with alcohol or betadine swabs.
- 4 Retract the underlying fascia will and drill a hole in the skull.

- 5 Insert a 26 gauge needle attached to a 2uL Hamilton syringe containing recombinant  $\alpha$ -syn PFF (5  $\mu$ g/2  $\mu$ L) or an equivalent volume of PBS into the striatum unilaterally (right hemisphere) at coordinates mediolateral 2.0 mm from bregma, anteroposterior -0.2 mm, and dorsoventral 2.8mm. at a rate of 0.4 ul per minute using a minisyringe pump (Model BS-800, Braintree Scientific).
- 6 Leave the needle in place for 5 minutes before withdrawing the needle.
- 7 Close the incision using sterile wound healing clips (7 mm stainless steel wound clip; Reflex 7 skin closure system) or sutures.
- 8 Immediately after the stereotaxic procedure, treat the stapled wound site with bupivacaine (0.4-0.8 ml/kg of a 0.25% solution, injected sc) or 0.5% xylocaine 1.4 ml/kg, injected sc, injected) or 0.5% xylocaine when necessary or as determined by animal use policy.