



DEC 15, 2022

## Scott's Tap Water Substitute

This protocol is published without a DOI.

COMMENTS 0

joshdm<sup>1</sup><sup>1</sup>UC Davis

joshdm

WORKS FOR ME

1

## ABSTRACT

Used as a bluing agent in H&amp;E staining

## PROTOCOL CITATION

joshdm 2022. Scott's Tap Water Substitute. **protocols.io**  
<https://protocols.io/view/scott-39-s-tap-water-substitute-ckhrut56>

## LICENSE

\_\_\_\_\_ This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

## CREATED

Dec 14, 2022

## LAST MODIFIED

Dec 15, 2022

## PROTOCOL INTEGER ID

74001

## PARENT PROTOCOLS

In steps of

[H&E staining](#)

## Measure and mix salts

1 weigh  20 g of  Magnesium Sulfate Heptahydrate **Fisher Scientific Catalog #M63-500**

2 weigh  2 g of  Sodium Bicarbonate **Sigma Aldrich Catalog #s-5761**

- 
- 
- 3 Mix in ddWater until completely dissolved  
store at room temperature

 8.3 to  8.7