



Jun 20, 2020

# Protocol for Exonuclease VIII, truncated (NEB #M0545)

New England Biolabs<sup>1</sup><sup>1</sup>New England Biolabs**1** Works for me This protocol is published without a DOI.**New England Biolabs (NEB)**Tech. support phone: +1(800)632-7799 email: [info@neb.com](mailto:info@neb.com) New England Biolabs Tech Support  
New England Biolabs

## ABSTRACT

Exonuclease VIII, truncated efficiently degrades linear dsDNA from 5' to 3' direction, leaving supercoiled dsDNA intact.

**\* Note:** For more precise results or partial digestions, we recommend titration of the enzyme to the intended substrate.

## EXTERNAL LINK

<https://neb.com/protocols/2019/07/02/protocol-for-exonuclease-viii-truncated-m0545>

## PROTOCOL CITATION

New England Biolabs 2020. Protocol for Exonuclease VIII, truncated (NEB #M0545). **protocols.io**  
<https://protocols.io/view/protocol-for-exonuclease-viii-truncated-neb-m0545-7sbhnan>

## EXTERNAL LINK

<https://neb.com/protocols/2019/07/02/protocol-for-exonuclease-viii-truncated-m0545>

## LICENSE

This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

## CREATED

Sep 27, 2019

## LAST MODIFIED

Jun 20, 2020

## OWNERSHIP HISTORY

Sep 27, 2019 Anita Broellochs protocols.io

Jun 18, 2020 New England Biolabs Tech Support New England Biolabs

## PROTOCOL INTEGER ID

28195

## MATERIALS

NAME	CATALOG #	VENDOR
EDTA	17892	Thermo Fisher
Exonuclease VIII truncated	M0545	New England Biolabs

## SAFETY WARNINGS

Please see SDS (Safety Data Sheet) for hazards and safety warnings.

1 

Set-up the reaction as follows:

COMPONENTS	50 µl REACTION
DNA	up to 1 µg
NEBuffer 4 (10X)	5 µl (1X)
Exonuclease VIII (truncated)	1 µl (10 units)
Nuclease-free H <sub>2</sub> O (NEB #B1500)	up to 50 µl

2 

Incubate at **37 °C** for **00:30:00**.

3 Stop reaction by adding *EDTA* to at least **11 Millimolar (mM)**.

4 Heat Inactivation **70 °C** for **00:30:00**.

5 To clean up treated samples, we recommend using **one of the following steps**:

- Column clean up (we recommend the Monarch<sup>®</sup>PCR & DNA Cleanup Kit, NEB #[T1030](#)) or
- Running the reaction on an agarose gel, and then extracting the DNA (we recommend the Monarch Gel Extraction Kit, NEB #[T1020](#)), or
- Performing a phenol/chloroform extraction followed by ethanol precipitation.