



Sep 06, 2020

u-Smell-it olfactory test

dkt¹¹D. Toomre, K Piironen

dkt : DRAFT - still under modification

In Development

This protocol is published without a DOI.



dkt

PROTOCOL CITATION

dkt 2020. u-Smell-it olfactory test . **protocols.io**
<https://protocols.io/view/u-smell-it-olfactory-test-bkxdxix6>



LICENSE

————— This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

Sep 06, 2020

LAST MODIFIED

Sep 06, 2020

PROTOCOL INTEGER ID

41669

GUIDELINES

1. Download u-Smell-it App from Apple Store or Android
2. Follow instructions

Alternative: Scan QR code on the Olfactory Card to Download

MATERIALS

NAME	CATALOG #	VENDOR
Multi-odorant olfactory test		

SAFETY WARNINGS

1. A positive result would need to be confirmed by a diagnostic test
2. A negative result may be inaccurate and miss people that have COVID

THIS IS NOT FDA approved and is only for research purposes.

DISCLAIMER:

DISCLAIMER – FOR INFORMATIONAL PURPOSES ONLY; USE AT YOUR OWN RISK

The protocol content here is for informational purposes only and does not constitute legal, medical, clinical, or safety advice, or otherwise; content added to protocols.io is not peer reviewed and may not have undergone a formal approval of any kind. Information presented in this protocol should not substitute for independent professional judgment, advice, diagnosis, or treatment. Any action you take or refrain from taking using or relying upon the information presented here is strictly at your own risk. You agree that neither the Company nor any of the authors, contributors, administrators, or anyone else associated with protocols.io, can be held responsible for your use of the information contained in or linked to this protocol or any of our Sites/Apps and Services.

BEFORE STARTING

Need:

1. Olfactory 'u-Smell-it' 5 window card
2. Need Cell Phone or Telephone

- 1 Receive Olfactory/Chemosensory test CARD
- 2 Scan QR code on the Olfactory Card (using any smart phone) and the camera function (they automatically recognize QR codes)

(Note: For people without access to a smart-phone there will be an automated call-in center that will provide voice prompts for choices of scent options and will connect to the centralized server, similar to the Apps)

- 3 Download 'u-Smell-it' App on iOS (Apple) or Android phone

- 4 Scan the QR-Code on the Card in the App

Note: After downloading the software Steps 1-3 are not necessary for additional tests (as would usually happen in repeated testing)

- 5 Follow the App instructions

5.1 Rub or scratch window #1 with finger(nail)

2s

5.2 In app choose the best of 4 scent options (one of which is the true scent) or pick the "no scent" option

- 6 Repeat 5.1 and 5.2 for all 5 scent window

- 7 Respond to symptom self-reporting questions

- 8 See the olfactory score (and once FDA approval is granted under an EUA will based on the test results will be advised to either: 1) self-isolate and seek medical advice (and follow up with a diagnostic test), 2) take another test (to improve the ability to resolve an ambiguous result such as a partial loss of smell) or 3) indicate that there the ability to smell is fine, but the subject should still exercise caution and that a negative result may miss positive people)

- 9 NEXT day, or