

Oct 01, 2021

Nuclei: Cytoplasm EKC Quantification Script (FIJI Macro)

Condon ND¹¹Institute for Molecular Bioscience, The University of Queensland, Brisbane, Australia

1 Works for me

Share

This protocol is published without a DOI.

Condon ND

ABSTRACT

This protocol describes how to download, install, and run the ImageJ/Fiji macro for quantifying the ratio fluorescence intensity of EKC in both the nuclei and cytoplasm. This protocol is used in the following eLife paper: <https://doi.org/10.7554/eLife.62196>

Briefly this macro takes directories of images and processes them in an autonomous batch-wise manner. Images should be 2-channel with channel one containing your EKC (or other tagged protein of interest) and channel 2 containing your nuclei signal. Nuclei are segmented and two selection regions are created (nuclei, and a surrounding band that is user definable). These two regions are measured for each identified nuclei and the mean fluorescence intensity is output into a single spreadsheet.

EXTERNAL LINK

<https://doi.org/10.7554/eLife.62196>

PROTOCOL CITATION

Condon ND 2021. Nuclei: Cytoplasm EKC Quantification Script (FIJI Macro). **protocols.io**
<https://protocols.io/view/nuclei-cytoplasm-ekc-quantification-script-fiji-ma-byp5pvq6>

FUNDERS ACKNOWLEDGEMENT

CZI
Grant ID: 2020-225648

MANUSCRIPT CITATION

 please remember to cite the following publication along with this protocol<https://doi.org/10.7554/eLife.62196>

KEYWORDS

Image Analysis, Quantification, ImageJ, FIJI, macro, EKC, signalling

LICENSE

This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

Oct 01, 2021

LAST MODIFIED

Oct 01, 2021

PROTOCOL INTEGER ID

53725

BEFORE STARTING

Having an up-to-date version of FIJI/ImageJ is critical for this macro. A fresh installation can be downloaded from [Fiji.sc](https://fiji.sc)

Downloading & installing the macro

- 1 Retrieve a copy of the macro file from either the [GitHub repo](#) or from the [Software Heritage archive](#)
 - 1.1 Download the .ZIP copy of the entire repo to your computer, and extract the compressed file.
 - 1.2 Launch FIJI/ImageJ on your computer.
 - 1.3 Install the plugin into your instance of FIJI by navigating to **Plugins > Macros > Install**
Choose to install the file called "Nuclei_cytoplasm-measuring_v1-5.ijm"

Running the macro

- 2 To launch the installed macro, navigate to **Plugins > Macros > Nuclei_cytoplasm-measuring_v1-5**
 - 2.1 The first window to appear when running the macro describes the Author, Licence and details the tasks the script will run.
 - 2.2 The second window warns the user that they will next be required to navigate to the directory which contains the image files. The user should also note the file extension of these files to enter into the upcoming preferences window.
 - 2.3 The preferences window prompts the user to input the File extension (for example .tif, .czi, .ism) as well as the banding distance (in pixels). The user can choose to run in batch mode (faster, background running) or not by selecting the relevant checkbox.
 - 2.4 The script will then run through each image within the chosen input directory, identifying nuclei, and measuring the EKC signal autonomously with no user input required.
 - 2.5 Upon completion the script will prompt the user with a dialog box stating "Put down that coffee! Your job is finished"

Reviewing the results

- 3 Navigate to the input directory location using your computers filesystem browser.
Select the newly created results directory called "Results<date&time>"
 - 3.1 For each image completed the following files are created:
<Filename>_cyto.zip [*this contains the ROIs for the cytoplasm measurements*]
<Filename>_green-selection.tif [*this image is the 32-bit thresholded image of Ch1*]
<Filename>_nuclei.zip [*this contains the ROIs for the nuclei measurements*]
<Filename>_nuclei-mask.tif [*this image is the thresholded mask of Ch2*]

<Filename>_selectionnumbers.png [*this image shows the ROI selection numbers corresponding to the quantification results*]


For the entire script run the following collated files are create:

Log.txt [*this is a copy of the Fiji Log window, with reported script version, settings used, directory location, and file completions*]

Results.xls [*results spreadsheet for quantification, see below*]

3.2 Results Spreadsheet information

Note if opening with Microsoft Excel you will see the following warning, this is because the file Results.xls was created outside of Microsoft Excel, click Yes to continue opening the document.



The following columns are provided in the Results.xls output file:

Filename: This is the filename for the input image.

Nuclei#: This is the nuclei ID (as per the corresponding output file
(<Filename>_selectionnumbers.png)

Nuclei Intensity: This is the measured Mean Intensity for the nuclei selection region

Cytoplasm Intensity: This is the measured Mean intensity for the cytoplasmic band selection region

Number of Nuclei: This is the total number of nuclei identified in this image (not it will repeat for each nuclei of this image)

Cyto/Nuclei: This is Cytoplasmic Intensity divided by the Nuclei Intensity

Nuclei/Cyto: This is the Nuclei Intensity divided by the Cytoplasmic intensity