



Jun 13, 2022

Worm Freezing Buffer



1

Adrien Assie¹, Buck Samuel¹¹Baylor College of Medicine

Adrien Assie: Uploader to Protocols.io, original protocol found on wormbook from Brenner et al. 1974

Buck Samuel: Uploader to Protocols.io, original protocol found on wormbook from Brenner et al. 1974

1

dx.doi.org/10.17504/protocols.io.14egnxjrzl5d/v1Adrien Assie
Baylor College of Medicine

Original protocol found on wormbook from Brenner et al. 1974

DISCLAIMER – FOR INFORMATIONAL PURPOSES ONLY; USE AT YOUR OWN RISK

The protocol content here is for informational purposes only and does not constitute legal, medical, clinical, or safety advice, or otherwise; content added to protocols.io is not peer reviewed and may not have undergone a formal approval of any kind. Information presented in this protocol should not substitute for independent professional judgment, advice, diagnosis, or treatment. Any action you take or refrain from taking using or relying upon the information presented here is strictly at your own risk. You agree that neither the Company nor any of the authors, contributors, administrators, or anyone else associated with protocols.io, can be held responsible for your use of the information contained in or linked to this protocol or any of our Sites/Apps and Services.

Worm Freezing buffer recipe according to Brenner, S. (1974). Genetics 77, 71. Found on wormbook see external link

DOI

dx.doi.org/10.17504/protocols.io.14egnxjrzl5d/v1http://www.wormbook.org/chapters/www_strainmaintain/strainmaintain.html#d0e920

Adrien Assie, Buck Samuel 2022. Worm Freezing Buffer. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.14egnxjrzl5d/v1>



protocol

Brenner, S. (1974). Genetics 77, 71



Samuel Lab Media and Buffers

protocol ,

Apr 10, 2019

Jun 13, 2022

22275

In steps of

Caenorhabditis elegans strain freezing procedure

Part of collection

Samuel Lab Media and Buffers

:




Original protocol found on wormbook from Brener et al. 1974

DISCLAIMER – FOR INFORMATIONAL PURPOSES ONLY; USE AT YOUR OWN RISK

The protocol content here is for informational purposes only and does not constitute legal, medical, clinical, or safety advice, or otherwise; content added to protocols.io is not peer reviewed and may not have undergone a formal approval of any kind. Information presented in this protocol should not substitute for independent professional judgment, advice, diagnosis, or treatment. Any action you take or refrain from taking using or relying upon the information presented here is strictly at your own risk. You agree that neither the Company nor any of the authors, contributors, administrators, or anyone else associated with protocols.io, can be held responsible for your use of the information contained in or linked to this protocol or any of our Sites/Apps and Services.

1 **150 mL** of **[M]100 % volume** glycerol

2 **25 mL** of **[M]1 Molarity (M)** Potassium Phosphate Buffer, pH 6.0

- 3  **2.85 g** of NaCl
- 4  **500 mL** of [M9 Buffer](#)
- 5 Adjust deionized water to  **1000 mL**
- 6 Autoclave