

JUN 01, 2023

## Modified Leeming and Notman Agar In 1 collection

Amy Gladfelter<sup>1</sup>

<sup>1</sup>Duke Universtiy



Amy Gladfelter

**ABSTRACT** 

Media recipe for growth of Malassezia species.

**ATTACHMENTS** 

Malassezia media\_1.pdf



## DOI:

dx.doi.org/10.17504/protocol s.io.kxygxm3ywl8j/v1

**Protocol Citation:** Amy Gladfelter 2023. Modified Leeming and Notman Agar. **protocols.io** 

https://dx.doi.org/10.17504/protocols.io.kxygxm3ywl8j/v1

License: This is an open access protocol distributed under the terms of the Creative Commons Attribution License, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working We use this protocol and it's working

Created: Nov 03, 2018

Last Modified: Jun 01, 2023

**PROTOCOL** integer ID:

17456

1 bacteriological peptone (Oxoid) 🚨 10 g

- 2 yeast extract A 2 g
- 3 desiccated ox bile (Oxoid) 4 8 g
- **4** glucose **△** 10 g
- 5 glycerol monostearate 🚨 0.5 g
- 6 agar 4 15 g
- 7 glycerol, 50% 🕹 20 mL
- 8 tween 60 A 5 mL
- 9 olive oil 4 20 mL

- 11 Sterilize by autoclaving.