



Version 1 ▼

Nov 04, 2022

OpenVent DNA Polymerase Production V.1

Jenny Molloy¹, Stephane Fadanka¹, Nadine Mowoh¹¹Beneficial Bio Ltd

1 Works for me

Share

This collection is published without a DOI.

[Reclone.org \(The Reagent Collaboration Network\)](#)[Beneficial Bio](#)

Jenny Molloy

ABSTRACT

This collection contains the protocols for production, quality control, packaging and use of Beneficial Bio's dehydrated OpenVent DNA Polymerase.

Format: 32 tubes (4x8-tube PCR strips) 20ul of Enzyme each

OpenVent DNA polymerase is an extremely high thermostable and high fidelity enzyme suitable for routine PCR applications and amplification of GC-rich or looped sequences. The enzyme is compatible with a wide range of templates and its robustness guarantees reliable amplification results in almost all PCR applications. The Eco and room temperature stable format - PCR reaction tube strips containing the dehydrated enzyme (10 reactions per tube), offers affordability and flexibility.

FEATURES

- Dehydrated formulation enables better stability at room temperature
- Ideal for a wide range of applications including GC-rich or looped sequences amplification
- Amplification of low copy DNA targets
- Eco formulation and flexible format

DOCUMENTS

[Product Manual](#)[Product Specification Sheet \(PSS\)](#)[Certificate of Analysis \(CoA\)](#)

COLLECTION CITATION

Jenny Molloy, Stephane Fadanka, Nadine Mowoh 2022. OpenVent DNA Polymerase Production. **protocols.io**
<https://protocols.io/view/openvent-dna-polymerase-production-ciucuesw>



LICENSE

————— This is an open access collection distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

IMAGE ATTRIBUTION

Beneficial Bio Ltd

CREATED

Nov 04, 2022

LAST MODIFIED

Nov 04, 2022

COLLECTION INTEGER ID

72292

ABSTRACT

This collection contains the protocols for production, quality control, packaging and use of Beneficial Bio's dehydrated OpenVent DNA Polymerase.

Format: 32 tubes (4x8-tube PCR strips) 20ul of Enzyme each

OpenVent DNA polymerase is an extremely high thermostable and high fidelity enzyme suitable for routine PCR applications and amplification of GC-rich or looped sequences. The enzyme is compatible with a wide range of templates and its robustness guarantees reliable amplification results in almost all PCR applications. The Eco and room temperature stable format - PCR reaction tube strips containing the dehydrated enzyme (10 reactions per tube), offers affordability and flexibility.



FEATURES



- Dehydrated formulation enables better stability at room temperature
- Ideal for a wide range of applications including GC-rich or looped sequences amplification
- Amplification of low copy DNA targets
- Eco formulation and flexible format



DOCUMENTS

[Product Manual](#)
[Product Specification Sheet \(PSS\)](#)
[Certificate of Analysis \(CoA\)](#)

FILES

- 

Production of cellular reagents using IPTG
Version 1
by Nadine Mowoh
- 

Functionality test (OpenVent polymerase, PCR Master Mixes)
Version 2
by Nadine Mowoh
- 

Nuclease Test (OpenVent polymerase, PCR Master Mix, DNA loading dye)
Version 1
by Nadine Mowoh