



VERSION 2

APR 05, 2024

🌐 Tail Suspension Test V.2

Marina Lorente Picón¹, Núria Peñuelas¹, Ariadna Laguna¹,
Marta Gonzalez-Sepulveda¹, Miquel Vila¹

¹Vall d'Hebron Research Institute

ASAP Collaborative Research Network

Vilalab Public



Miquel Vila
VHIR-CIBERNED-ASAP

ABSTRACT

Tail Suspension Test for mice

OPEN  ACCESS



DOI:

dx.doi.org/10.17504/protocols.io.6qpvr3222vmk/v2

Protocol Citation: Marina Lorente Picón, Núria Peñuelas, Ariadna Laguna, Marta Gonzalez-Sepulveda, Miquel Vila 2024. Tail Suspension Test. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.6qpvr3222vmk/v2> Version created by [Miquel Vila](#)

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: Apr 05, 2024

Last Modified: Apr 05, 2024

PROTOCOL integer ID: 97833

Funders Acknowledgement:

Aligning Science Across
Parkinson's

Grant ID: ASAP-020505

- 1 Suspend animals by their tails with tape in a suspension bar.

- 2 To avoid the tail climbing behavior, pass a 2cm methacrylate tube through the tail before suspending the animal.

- 3 Quantify the escape-oriented behaviors (i.e. fore and hind limbs movement) during six minutes. Vocalizations were registered during this period with a yes/no score.

- 4 Calculate the total immobilization time (s) as the sum of all the time the animal was not performing any escape-oriented behavior.