




SEP 28, 2023

## Sea Water - Blue Treasure

Willian Barela

Costa<sup>1,2</sup>

<sup>1</sup>FabLab Casa Firjan; <sup>2</sup>@casafirjan



Willian Barela Costa

### ABSTRACT

Sea water production for microalgae cultivation using the Blue Treasure supplier.

25 - 35 Salinity

[M] 27 mg/mL

### MATERIALS

Salt from Blue Treasure product

Distilled water

OPEN ACCESS



DOI:

[dx.doi.org/10.17504/protocols.io.36wgq36b5lk5/v1](https://dx.doi.org/10.17504/protocols.io.36wgq36b5lk5/v1)

**Protocol Citation:** Willian Barela Costa 2023. Sea Water - Blue Treasure. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.36wgq36b5lk5/v1>




**License:** This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working  
We use this protocol and it's working

**Created:** Sep 26, 2023

**Last Modified:** Sep 28, 2023

## Salting process - 1 liter

- 1  0.027 g/mL ->  Sample 27 g Sea Sal to  1 L Water distilled
- 2 Shake until diluted
- 3 Make sure the salinity is within the desired range using the salinity refractometer. In this protocol, the focus is 27 salinity.