



VERSION 2

FEB 10, 2023

OPEN ACCESS

DOI:

[dx.doi.org/10.17504/protocols.io.eq2lynymvx9/v2](https://dx.doi.org/10.17504/protocols.io.eq2lynymvx9/v2)

**Protocol Citation:** Malcolm A Barnard, Stephen Powers 2023. Disposal of Environmental Samples from Water Bodies Infested with Zebra Mussels . **protocols.io** <https://dx.doi.org/10.17504/protocols.io.eq2lynymvx9/v2> Version created by **Malcolm A Barnard**

**License:** This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working  
We use this protocol and it's working

**Created:** Aug 20, 2022

**Last Modified:** Feb 10, 2023


**PROTOCOL integer ID:**  
68955

## 🌐 Disposal of Environmental Samples from Water Bodies Infested with Zebra Mussels V.2

Malcolm A Barnard<sup>1</sup>, Stephen Powers<sup>1</sup>

<sup>1</sup>Baylor University

Powers Lab Baylor



Malcolm A Barnard

Baylor University

### DISCLAIMER

DISCLAIMER – FOR INFORMATIONAL PURPOSES ONLY; USE AT YOUR OWN RISK

The protocol content here is for informational purposes only and does not constitute legal, medical, clinical, or safety advice, or otherwise; content added to [protocols.io](#) is not peer reviewed and may not have undergone a formal approval of any kind. Information presented in this protocol should not substitute for independent professional judgment, advice, diagnosis, or treatment. Any action you take or refrain from taking using or relying upon the information presented here is strictly at your own risk. You agree that neither the Company nor any of the authors, contributors, administrators, or anyone else associated with [protocols.io](#), can be held responsible for your use of the information contained in or linked to this protocol or any of our Sites/Apps and Services.

### ABSTRACT

Zebra Mussels (*Dreissena polymorpha*) are a highly invasive species that have been invading water bodies around the United States and beyond. This protocol for sample disposal for sample collected from know zebra mussel infested water bodies to prevent spread though local water treatment facilities by bleaching the samples to kill any larvae or other potential zebra mussels in the sample water.

### MATERIALS

Bleach (10% Sodium Hypochlorite)  
Dishwashing bin to add sample water and bleach

## SAFETY WARNINGS



Do not mix acids and bleach. This will create chlorine gas.

Wear gloves, eye protection, and a lab coat while performing this protocol.

## Sample Killing

4h 5m

1 Pour all environmental sample water from bottles into a dish pan.

2 Save the bottles for acid washing

### Safety information

Do not add bleach to bottles being cleaned using an acid bath and acid as the combination of acid and bleach produces chlorine gas.

3 Add 50 mL of bleach per 1 L of sample water.

4 Let the bleached sample sit for a minimum of four hours if not until the next day.

4h

## Sample Disposal

5m

5 Dispose of bleached sample down drain with excess water or other waste pathway as described by local environmental health and safety guidelines.

6 Save bin for future sample disposal or bleach work.

