



Jul 21, 2021

♠ An Optimized Protocol for the Generation of Midbrain Dopamine Neurons from human iPSCs under Defined Conditions


Carlos W. Gantner^{1,2,3,4,5}, Agustín Cota-Coronado^{1,6,2}, Lachlan H. Thompson¹, Clare L Parish^{1,7,8}

¹The Florey Institute of Neuroscience and Mental Health, The University of Melbourne, Parkville, VIC 3052, Australia;

³Present address: Mammalian Embryo and Stem Cell Group, Department of Physiology, Development and Neuroscience, University of Cambridge, Cambridge CB2 3DY, UK;


⁶Biotecnología Médica y Farmacéutica CONACYT, Centro de Investigación y Asistencia en Tecnología y Diseño del Estado de Jalisco (CIATEJ), Guadalajara 44270, Mexico;

⁷Lead Contact; ⁸Correspondence: clare.parish@florey.edu.au



²These authors contributed equally;

⁴Technical Contact; ⁵Correspondence: cg731@cam.ac.uk;



Here, we describe a xeno-free, feeder-free, and chemically defined protocol for the generation of ventral midbrain dopaminergic (vmDA) progenitors from human pluripotent stem cells (hPSCs). This simple-to-follow protocol results in high yields of cryopreservable dopamine neurons across multiple hPSC lines. Wnt signaling is the critical component of the differentiation and can be finely adjusted in a line-dependent manner to enhance production of dopamine neurons for the purposes of transplantation, studying development and homeostasis, disease modeling, drug discovery, and drug development.

For complete details on the use and execution of this protocol, please refer to Gantner et al. (2020) and Niclis et al. (2017a).

HIGHLIGHTS

- Reproducible differentiation of human dopamine neurons from multiple hPSC lines
- Dopamine progenitors can be cryopreserved for downstream applications
- Dopamine neurons mature in vitro, enabling screening or developmental studies
- Transplanted dopamine progenitors are capable of restoring motor function

ATTACHMENTS

An Ontimized Protocol for the Generation of Midbrain Dopamine Neurons under Defined Conditions.pdf

DOI

dx.doi.org/10.17504/protocols.io.bu7ynzpw

PROTOCOL CITATION

Carlos W. Gantner, Agustín Cota-Coronado, Lachlan H. Thompson, Clare L Parish 2021. An Optimized Protocol for the Generation of Midbrain Dopamine Neurons from human iPSCs under Defined Conditions. protocols.io

https://dx.doi.org/10.17504/protocols.io.bu7ynzpw

A

m protocols.io

07/21/2021

2 Citation: Carlos W. Gantner, AgustÃân Cota-Coronado, Lachlan H. Thompson, Clare L Parish (07/21/2021). An Optimized Protocol for the Generation ofÃâ Midbrain Dopamine Neurons from human iPSCs under DefinedÃÂ Conditions. https://dx.doi.org/10.17504/protocols.io.bu7ynzpw

Gantner et al., An Optimized Protocol for the Generation of Midbrain Dopamine Neurons under Defined Conditions, STAR Protocols (2020), https://doi.org/10.1016/j.xpro.2020.100065

LICENSE

This is an open access protocol distributed under the terms of the Creative Commons Attribution License, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

May 21, 2021

LAST MODIFIED

Jul 21, 2021

OWNERSHIP HISTORY



PROTOCOL INTEGER ID

50136

GUIDELINES




Figure 1. Differentiation Overview and Morphological Changes during vmDA Differentiation

(A) Schematic representation of the differentiation protocol indicating when each factor is added to the base media to appropriately pattern hPSC toward the vmDA lineage.

(B) Morphological changes during differentiation. At day 3, the cells are a flat, uniform monolayer. By day 11, the density has increased, and the wells appear darker under the microscope. By days 18–20, fibers begin to appear as neurons accumulate and aggregate together and by day 25 vmDA neurons display extended axonal processes. Scale bar, 100 µm.




Figure 2. Characterization of Developing vmDA Progenitors and Neurons with 2.5 mM CHIR99021
(A) Immunofluorescence analysis of developing vmDA progenitors differentiated from either H9:PITX3-eGFP or RM3.5:PITX3-eGFP using 2.5 mM CHIR99021. At day 13, vmDA progenitors should homogenously express OTX2, FOXA2 and Nestin, indicating correct midbrain floorplate lineage acquisition.

(B) At day 25, vmDA neurons expressing FOXA2, PITX3, and TH should be readily observed if differentiation has been successful. Scale bar, 100 mm.

A	В	С	D	Е	F	G	Н	I	J	K
	Nestin	OTX2	FOXA2	LMX1A	EN1	NURR1	PITX2	BARHL1	(PITX3)eGFP	TH
	1:1000	1:1000	1:500	1:200	1:50	1:200	1:500	1:500	1:1000	1:1000
Day 13	+++	+++	+++	++	+/-	-/+	+/-	+/-	-	-
Day 16	+++	+++	+++	++	++	++	+/-	+/-	+/-	+/-
Day 25	++	+++	+++	++	+/-	+	+/-	+/-	++	++
Day 60	+/-	++	++	+	+/-	+/-	+/-	+/-	+++	+++

Table 5. Expression of vmDA Genes during Differentiation Expression level: - none; +/- low; + moderate; ++ high; +++ very high.




Figure 3. vmDA Progenitors Can Be Cryopreserved and Thawed with High Viability

(A) Differentiation schematic highlighting two time points where vmDA progenitors can be frozen and/or transplanted with high viability.

(B) Live/dead analysis of DNA content using FACS highlights that cells are viable after cryopreservation at days 13 or 19.
(C) By day 25, thawed and replated vmDA progenitors mature into FOXA2+/PITX3-eGFP+ vmDA neurons in vitro. Scale bar, 100

EXPECTED OUTCOMES

One day after seeding, undifferentiated hPSC should tessellate the culture surface such that there are no (or very small) gaps between cells. As differentiation proceeds cells will actively proliferate and appear denser, until individual cells can no longer be resolved under the microscope (typically day 5). At the early stages of differentiation there should be very little cell death. By day 13 the majority of the culture should appear homogenous in distribution and density across the culture well, co-express FOXA2 and OTX2 (≥80%) and will begin to appear darker under the microscope as cell density increases (Figures 1B and 2; Table 5). Note, some death may occur following the transition from Induction to Maturation Media. Addition of the g-secretase inhibitor DAPT, which drives cell cycle exit, often leads to clear neuronal processes by day 20, although this is cell line-dependent and does not necessarily reflect the final TH+ vmDA neuron yield.

By day 20 TH+ (or PITX3-eGFP+) neurons are readily observable and are increasingly observed until peaking at

day 40 (D20: \sim 5%; D25: \sim 25%; D40: \sim 50%) (Figure 2; Table 5; See Niclis et al. (2017a) for more detailed quantification across numerous hPSC lines). vmDA neurons should also co-express FOXA2, OTX2, NURR1, EN1, LMX1A (Figure 2; Table 5). vmDA progenitors and neurons will increasingly clump together after day 30 as they continue to mature. At day 40 vmDA neurons in culture are capable of synthesizing and metabolizing dopamine, and by day 60 vmDA neurons display characteristic vmDA electrophysiological properties (Niclis et al., 2017a).

For the purpose of transplantation, we typically obtain total cell yields of $1-2 \times 10^6/\text{cm}^2$ at time points amenable to implantation (days 13-22). Demonstrative of the functionality and integration of these cells, we observe correction in the unilateral 6-hydroxydopamine rat model of DA depletion/Parkinson's disease in amphetamine-induced rotational asymmetry after 12-16 weeks.

LIMITATIONS

We have successfully utilized this protocol to generate vmDA neurons from multiple hESC and hiPSC cell lines (including H9, HES3, RUES2, 409B2 parental lines and several subclones). However, differentiation efficiency does vary between hPSC lines. We have found that the majority of this variation is Wnt-dependent. As such, we recommend optimization of the CHIR99021 concentration for each cell line. Typically, titration between $2-3~\mu M$ is sufficient to optimize vmDA yield (Figure 4). This is particularly important for iPSC-based disease modeling studies, that typically derive numerous cell lines. Even under optimal conditions, where OTX2+/F0XA2+, putative vmDA progenitors make up \geq 90% of the culture, only \sim 50% of these will mature into TH+ vmDA neurons. This is likely due to the aberrant maintenance of an immature state in vitro or the generation of adjacent midbrain lineages. Typically, we do not observe significant 5-HT, GABA or Glutamate (\geq 2%) contamination within vmDA cultures. If specific analysis of the vmDA population is desired (e.g. for in vitro disease modeling of vmDA phenotypes, vmDA-centric drug development, or grafting in vivo), we suggest the use of an appropriate fluorescent reporter (such as PITX3-GFP or TH-GFP) to allow for isolation (e.g., for multi-omic analyses) and/or tracking in vitro or after transplantation. See Niclis et al. (2017b) for an example of how a fluorescent reporter can augment analysis.

Finally, the use of fully defined media, utilizing xenogeneic-free reagents (amenable to clinical translation), results in less variability and better efficiency between differentiations and cell lines (Niclis et al., 2017a). However, the use of animal product-free reagents can add unnecessary expense when such stringency may not be required. For other applications we recommend using the more traditional media formulations, which are typically less expensive and work equally well. We regularly interchange xenogeneic and xenogeneic-free reagents (e.g. mTeSR2/mTeSR1) without altering the protocol, although this should be determined on a case-by-case basis by the end user




Figure 4. CHIR99021 Concentration Determines the Rostrocaudal Positioning of Developing Neural Progenitors

Titrating CHIR99021 is advisable for each individual cell line.

(A) Neural progenitors at day 13 exposed to increasing levels of [CHIR99021] (0.5–3.0 μ M) demonstrating a loss of fore/midbrain marker OTX2 and increased EN1 at higher concentrations. PITX2 and BARHL1 are expressed in the rostral midbrain and are helpful to fine-tune CHIR99021 concentration. We have found that 2.5 mM consistently generated vmDA progenitors from multiple hPSC cell lines.

(B) PITX3-eGFP+/TH+ vmDA neurons are born only after appropriate CHIR99021 addition and are progressively lost when [CHIR99021] is too high and cultures are overcaudalized. Scale bar, 100 μm

(iii) protocols.io 5 07/21/2021

FOR I ROUBLESHOUTING, RESOURCE AVAILABILITY, ACKNOWLEDGMENTS, AUTHOR CONTRIBUTIONS, DECLARATION OF INTERESTS and REFERENCES, please take a look at the attachment.

MATERIALS TEXT

KEY RESOURCES TABLE

Α	В	С
REAGENT or RESOURCE	SOURCE	IDENTIFIER
Antibodies	!	
Rabbit polyclonal anti-	Novus Biologicals	NBP1-86513,
BARHLI1		RRID:AB_11034569
Mouse monoclonal anti- EN1	DSHB	4G11, RRID:AB_528219
Goat polyclonal anti-FOXA2	Santa Cruz	sc-6554,
	Biotechnology	RRID:AB_2262810
Chicken polyclonal anti-GFP	Abcam	ab13970,
		RRID:AB_300798
Rabbit polyclonal anti-	Millipore	AB10533,
LMX1A		RRID:AB_10805970
Rabbit polyclonal anti-	Santa Cruz	Sc-990, RRID:AB_2298676
NURR1/NUR77	Biotechnology	
Goat polyclonal anti-OTX2	R&D Systems	AF1979,
		RRID:AB_2157172
Sheep polyclonal anti-	R&D Systems	AF7388,
PITX2		RRID:AB_11128639
Rabbit polyclonal anti-TH	Pel-Freez Biologicals	P40101-0,
		RRID:AB_461064
Sheep polyclonal anti-TH	Pel-Freez Biologicals	P60101-0,
Observiceds Doubids and Do		RRID:AB_461070
Chemicals, Peptides, and Re		AT404
Accutase®	Innovative Cell	AT104
A	Technologies	A 4 4 0 0
Ascorbic acid	Sigma-Aldrich	A4403
B27 TM Supplement or	Gibco	17504044
CTS TM B27 TM Supplement	0.1	A1486701
B27 TM Supplement -vitA or CTS TM B27	Gibco	12587010
Supplement -vitA		A3353501
Brain-derived neurotrophic	R&D Systems	248-BDB
factor	Rad Systems	240-000
CHIR99021	Miltenyi Biotech	130-103-926
Dibutyryl cAMP	Tocris Bioscience	1141
DMEM/F-12 or CTS TM	Gibco	11320033
DMEM/F-12 01 C13****	GIDCO	A1370801
Dimethyl sulfoxide	Sigma-Aldrich	D2650
Glial cell line-derived	-	
neurotrophic factor	R&D Systems	212-GD
<u>'</u>	O:h	25050061
GlutaMax TM Supplement	Gibco	35050061

Insulin-Transferrin-	Gibco	51300044
Selenium-Sodium		
Pyruvate (ITS-A)	0.1	10000010
KnockOut TM serum replacement or CTS TM	Gibco	10828010 12618012
KnockOut TM serum		12010012
replacement		
Laminin-521	BioLamina	LN521
Laminin - mouse	Sigma-Aldrich	L2020
LDN193189	Miltenyi Biotech	130-103-925
MEM Non-Essential Amino	Gibco	11140035
Acids		
mTeSR TM 1 or TeSR TM 2	STEMCELL Technologies	85880
	-	05860
N2 Supplement or CTS TM	Gibco	17502048
N2 Supplement		A1370701
N-[N-(3,5-	Sigma-Aldrich	D5942
Difluorophenacetyl)-L-		
alanyl]-S-phenylglycine t-		
butyl ester (DAPT)		
Neurobasal TM or CTS TM	Gibco	A3582901
Neurobasal TM		A1371201
PBS ^{+/+} (with	Gibco	13492609
Calcium/Magnesium)		
PBS ^{-/-}	Gibco	14190144
Pierce TM 16%	Thermo Fisher	28906
formaldehyde		
Purmorphamine	Miltenyi Biotech	130-104-465
ReLeSR TM	STEMCELL Technologies	05873
SB431542	R&D Systems	1614
SHH C25II	R&D Systems	464-SH
TGFb3	Peprotech	100-36E
Trypan blue	Sigma-Aldrich	T8154
Y27632 (ROCK inhibitor) or	Tocris Bioscience	1254
Y27632-GMP		TB1254-GMP
Critical Commercial Assays		
MycoAlert TM Mycoplasma	Lonza	LT07-218
Detection Kit		
Experimental Models: Cell Lin	nes	
H9 PITX3-eGFP hESC	Colin Pouton Laboratory.	N/A
	(Watmuff et al., 2015)	
RM3.5 PITX3-eGFP hiPSC	In house	N/A
Other		
48-well plate - Nunclon	Thermo Fisher	150687
Delta surface		
Bright-Line TM	Sigma-Aldrich	Z359629
hemocytometer		
Cryovials – 1.5 mL	Thermo Fisher	366656
Mr. Frosty TM freezing	Thermo Fisher	5100-0001
container		
T25 Flask	Corning	353108



MATERIALS AND EQUIPMENT

Laminin-521 Dilution (Tables 1 and 2)

Α	В	С	D	
	Laminin-521	PBS+/+ (μL)	Final Volume	
	(μL)		(μL)	
T25	150	2850	3000	
48-well plate	7.5	142.5	150	
96-well plate	3.5	66.5	70	

Table 1. Laminin-521 Dilution for 5 μg/mL

Α	В	С	D
	Laminin-521	PBS+/+ (μL)	Final Volume
	(μL)		(μL)
T25	300	2700	3000
48-well plate	15	135	150
96-well plate	7	63	70

Table 1. Laminin-521 Dilution for 10 μg/mL

Base Media Composition (Tables 3 and 4)

Α	В
Induction media (N2B27)	Volume (mL)
DMEM/F12	47.75
Neurobasal media (NBM)	47.75
B27 -Vitamin A (1x)	2
N2 (1x)	1
Glutamax (1x)	1
Penicillin/Streptomycin (0.5x)	0.5
Total	100 mL

Table 3. Induction Media Composition

A	В
Maturation Media (NBB27)	Volume (mL)
DMEM/F12	46.75
Neurobasal media (NBM)	46.75
B27 +Vitamin A (1x)	2
N2 (1x)	1
Non-essential amino acids (1x)	1
ITS-A (1x)	1
Glutamax (1x)	1
Penicillin/Streptomycin (0.5x)	0.5
Total	100 mL

Table 4. Maturation Media Composition

SAFETY WARNINGS

For hazard information and safety warnings, please refer to the SDS (Safety Data Sheet).

Note: All steps should be carried out under sterile conditions and according to local regulatory requirements (e.g. Australian Std 2243 requiring all human tissues to be confined to Class II biosafety cabinets)

Note: Cells are cultured throughout at § 37 °C in [M]5 % CO2.

BEFORE YOU BEGIN: Pluripotent Stem Cell Thawing

2h 18m

1 🔲 /

2h

Coat a T25 flask with [M]5 μ g/mL - [M]10 μ g/mL Laminin-521 and leave in the incubator for $\geq \bigcirc$ 02:00:00 (see Tables 1 and 2 for laminin-521 dilution).

Note: [M]5 μ g/mL of Laminin-521 is sufficient for cell lines that have been adapted to Laminin-521. For non-adapted lines (i.e. hPSC lines maintained on substrates other than laminin-521) [M]10 μ g/mL may be required to ensure attachment during initial passages. Laminin-521 can also be added for 12:00:00 - 18:00:00 at 4 °C, although 20:00:00 at 37 °C is preferable.

CRITICAL: Laminin-521 must be diluted in PBS+/+. Diluting in PBS^{-/-} will result in poor attachment and survival. Ensure that enough laminin solution is made to coat the entire flask throughout the coating period. Do not use flasks which have dried as cells will not attach in these areas.

2

15m

Equilibrate \blacksquare 15 mL pre-made mTeSR1 in a tightly closed 15 mL tube in the incubator for $\ge \bigcirc 00:15:00$.

3

Rapidly thaw a cryovial of hPSC in a § 37 °C water bath and transfer the contents to

■10 mL pre-warmed mTeSR1 solution using a P1000 pipette.

4 🚇 🥂

3m

Centrifuge cells at 300 x g, 21°C, 00:03:00 , aspirate supernatant, and resuspend the cell pellet in

■5 mL mTeSR1 media, supplemented with [M]10 Micromolar (μM) Y27632.

CRITICAL: Addition of the ROCK inhibitor Y27632 is essential for hPSC survival, particularly if the cells have been cryopreserved as single cells rather than aggregates.

protocols.io

07/21/2021

Aspirate laminin from the T25 flask and immediately plate the cell suspension, gently rocking the flask multiple times up-and-down and then side-to-side to evenly distribute the cells.

BEFORE YOU BEGIN: Pluripotent Stem Cell Culture



One day after thawing, wash once with PBS^{-/-} to remove debris, and replace with fresh media.

7 Continue to culture hPSC until they reach 70%–80% confluence, changing media every day.

Note: Generally, cells will recover quickly after thawing and should reach confluency within 5 days. However, if cells are growing slowly it is best to passage at a lower density (50% - 60%) confluence or after 5-6 days.

8 8

15m

To passage hPSC, wash once in PBS^{-/-} and add **□2 mL ReLeSR**. Leave at § **Room temperature** and monitor cell morphology under the microscope until colonies are loosely attached (© 00:05:00 − © 00:10:00).

9 Remove ReLeSR carefully and rapidly tap the flask hard 3-4 times on each long edge to detach colonies.

10

Resuspend by dissociating colonies in **2 mL equilibrated mTeSR1** using a P1000 pipette. Dilute cells (1:10–1:20) and transfer using a 5 mL serological pipette to a Laminin-521 coated T25 flask.

Note: Morphology during passage can vary between cell lines and the optimal timing of ReLeSR should be determined accordingly. The majority of cells should remain as 10–50 cell aggregates, removing the need for Y27632. Add Y27632 if passaged as single cells. Likewise, the growth rate of each line will differ, and passage ratios should be adjusted to compensate. We recommend passaging every 3–5 days.

Optional: Cryopreserve stem cells at low passage number by adding [M]10 % DMSO, [M]45 % KnockOut Serum Replacement and [M]45 % mTeSR to isolated cell pellets (after ReLeSR dissociation) and slow-freeze at & -80 °C . 3-4 cryovials can be generated from a single, confluent T25 flask. Transfer to liquid nitrogen for long-term storage.

Pause Point: hPSC can be maintained by passaging every 3–5 days until the desired experimental start time. Thinking ahead and knowing the growth characteristics of the cell line are helpful in aligning cells at 70%–80% confluence with the desired start date.

CRITICAL: Always passage thawed cells ≥ 2 times before beginning a differentiation experiment, or until the cell line displays normal morphology and growth rate characteristics.

CRITICAL: Routine mycoplasma testing is essential to maintaining high-quality hPSC stocks that will consistently differentiate. The quality of stem cells will dictate differentiation efficiency and should be closely monitored by testing, morphology and immunostaining (e.g. \geq 95% OCT4+). In addition, we do not recommend using high passage number hPSC.

BEFORE YOU BEGIN: Preparing Differentiation Media

20m

11 Prepare media as outlined in the Materials and Equipment section below (Tables 3 and 4).

Note: Always prepare media as close to use as possible and use within 2 weeks. It is best to make Induction Media the day of the experiment and Maturation Media at day 11. Where possible, add each factor (SHH, CHIR, etc.) to the base media on the day to minimize degradation.

Plating hPSC for Differentiation

1h

12 Timing: 1 h

This section describes the passaging of hPSC cultures into monolayers appropriate for differentiation.

13



2h

Coat a 48-well plate with [M]5 μ g/mL Laminin-521 and leave in the incubator for $\geq \bigcirc$ 02:00:00 (see Tables 1 and 2 for laminin-521 dilution).

Note: Differentiation can proceed in any size multi-well plate or flask. Adjust the Laminin-521 dilution, seeding density, and media volume accordingly.

14



10m

Passage hPSC at 70%–80% confluency by washing once in PBS^{-/-}, adding $\blacksquare 2$ mL Accutase, and placing in the incubator for 00:04:00 – 00:06:00. Monitor the morphology of the cells closely during this period. Once colonies show a characteristic bright halo and are close to detaching, gently pipette using a P1000 until $\ge 90\%$ of the flask has detached.

15



3m

Rapidly dilute Accutase solution in ≥ □3 mL mTeSR1 and centrifuge at ⊗300 x g, 21°C, 00:03:00.

16



Aspirate supernatant, resuspend the cell pellet in fresh § 37 °C warmed mTeSR1 supplemented with [M] 10 Micromolar (µM) Y27632 and count live cells using a hemocytometer.

Note: Resuspend hPSC at 1-2 x 106 cells/mL to enable rapid, accurate counting. Typically, we resuspend hPSC in 1 mL/5-10 cm². Each cell line grows at a different density and therefore the correct volume to resuspend cells in should be determined individually.



25m

Seed $0.3 \times 10^6 \text{ cells/cm}^2$ (48 well surface area = $\sim 1 \text{ cm}^2$) in $\geq 400 \text{ µl mTeSR1} + \text{Y27632}$. Rock the plate in a "figure eight" movement to evenly distribute cells and leave to settle for © 00:10:00 - © 00:15:00 before placing in the incubator.

Note: It is crucial to start the differentiation with an even monolayer of cells. Carefully dissociating to a single cell solution is important as it allows for an accurate count and an even distribution upon seeding. It is equally important not to prolong exposure of the hPSCs to enzymatic dissociation as this will impair survival. The duration of enzymatic treatment and the balance between optimal dissociation and survival may vary between cell lines.

Note: Initial seeding density should be optimized for each cell line to ensure confluency 24 h after seeding. Typically, we seed 0.2-0.4 x 10⁶ cells/cm². Trypan blue exclusion dye can be utilized during cell counting to ensure that only viable cells are included in the final cell seeding calculation.

CRITICAL: Seeding too many cells can lead to early detachment during differentiation, while seeding an insufficient number of cells can lead to poor survival and off-target differentiation.

vmDA Progenitor Differentiation

Timing: 13 days 18

This step describes the steps needed to differentiate hPSCs into vmDA progenitors.

See Figure 1A for a schematic of the differentiation.



Day 0: 24 h after seeding check under the microscope that hPSC have survived and appear as a single, confluent (\geq 95%) monolayer with no obvious differentiation or uneven distribution.

Note: No (or very few small) gaps should be present the day after seeding the cells. If it is hard to see, wash twice with PBS^{-/-} and check under the microscope. In addition, there should be very little death. If there is a large amount of death or cells are not confluent, it is best to discard the cells and start again.

- Prepare Induction Media and supplement with [M]10 Micromolar (μM) SB431542 and [M]200 Nanomolar (nM) LDN193189, SMAD inhibitors of the Transforming Growth Factor-b and Bone Morphogenetic Protein pathways, respectively.
- 21

Gently wash each well twice with PBS^{-/-} to remove residual mTeSR and replace media with $\blacksquare 500 \ \mu L$ supplemented Induction Media.

- Day 1: Change media to 500 μL Induction Media supplemented with [M]10 Micromolar (μM) SB431542, [M]200 Nanomolar (nM) LDN193189, as well as ventralizing factors [M]100 ng/mL C25II SHH and [M]2 Micromolar (μM) purmorphamin.
- Days 2, 3 and 4: Change media to 500 μl supplemented with [M]10 Micromolar (μM) SB431542, [M]200 Nanomolar (nM) LDN193189, [M]100 ng/mL C25II SHH, [M]2 Micromolar (μM) purmorphamine and [M]2.5 Micromolar (μM) Wnt-agonist CHIR99021 to caudalize developing progenitors.

Note: Media is changed every 24 h and 500 mL/48-well is added at days 0-4. At day 5 and beyond media is changed every other day and the volume is increased to 1 mL/48-well.

- Day 5: Change media to 1 mL Induction Media supplemented with [M]200 Nanomolar (nM) LDN193189, [M]100 ng/mL C25II SHH, [M]2 Micromolar (μM) purmorphamine and [M]2.5 Micromolar (μM) CHIR99021.
- Days 7 and 9: Change media to 1 mL Induction Media supplemented with

 [M]200 Nanomolar (nM) LDN193189 and [M]2.5 Micromolar (μM) CHIR99021.
- 26 ₺

Day 11: Wash cultures twice in PBS^{-/-} and switch base media to Maturation Media with [M]20 ng/mL rhBDNF, [M]20 ng/mL rhGDNF, [M]0.1 Milimolar (mM) dibutyryl cAMP, [M]200 Nanomolar (nM) Ascorbic Acid, [M]1 ng/mL rhTGFb3, and [M]10 Micromolar (μM) DAPT (Maturation Media + BGCATD), supplemented with [M]2.5 Micromolar (μM) CHIR99021.

Optional: Analyze cells at the end of progenitor patterning, day 13, for markers of vmDA progenitors by fixation and immunofluorescence or collection of cells for RT-qPCR. See below for the immunostaining protocol. By day 13 the majority (≥80%) of cells should co-express OTX2, FOXA2 and Nestin (<u>Figure 2A</u>; <u>Table 5</u>). Together, these markers indicate correct specification of ventral midbrain progenitors.

Note: If vmDA progenitors are too dense at day 13 for accurate quantification, it is possible to passage the cells with Accutase. If so, replate cells in Maturation Media + BGCATD, supplemented with [M] **10 Micromolar (\muM) Y27632** at 0.3–0.45 x 10⁶/cm². Fix 24 h later, after cells have attached, for analysis or continued culture.

CRITICAL: Developing neural progenitors are extremely sensitive to CHIR99021 concentration. Take care when diluting new vials of CHIR99021 and precisely pipette the correct volume during the experiment. It is best to always add CHIR99021 on the day if possible.

Terminal Differentiation and Long-Term Culture of vmDA Neurons

27 Timing: 12+ days

This step describes the long-term culture of vmDA progenitors and the generation of vmDA neurons.

28



Days 13-30. Change media every other day to Maturation Media + BGCATD.

Optional: Transplantation of vmDA progenitors is optimal between days 13–22. Wash cells twice in PBS^{-/-}, add 400 mL/48-well of Accutase and incubate for ⊙ 00:07:00 – ⊙ 00:10:00 at § 37 °C. Monitor the wells closely under a microscope and test detachment every few minutes by gently pipetting the Accutase directly onto the cells. Once cells begin to detach, gently pipette using a P1000 to collect ≥90% of the well and further dissociate vmDA progenitors. Accutase can be left on for up to 20 min for slowly detaching, more dense cultures. Arrest enzymatic digestion when cells have dissociated into small (5–20 cell) aggregates by adding Maturation Media + BGCATD, supplemented with [M]10 Micromolar (µM) Y27632 and centrifuge at

300 x g, 21°C, 00:03:00 . Aspirate the supernatant, resuspend the cell pellet in (1 mL/cm²) Maturation Media + BGCATD + [M]10 Micromolar (μM) Y27632 and count viable progenitors using a hemocytometer and Trypan blue exclusion. Centrifuge at 300 rpm, 21°C, 00:03:00 , aspirate supernatant, and resuspend at 0.1–0.2 x 10⁶/μL in Maturation Media + BGCATD + [M]10 Micromolar (μM) Y27632 and store 8 On ice prior to transplantation.

Note: From day 18–20, vmDA neurons will become increasingly abundant and can be indirectly observed by the formation of neuronal bunches and axonal projections. By day 25, these neuronal bundles become dense and highly arborized (**Figure 1B**), making dissociation challenging and resulting in poor survival after plating and following transplantation. Viability when counting cells for further culture or transplantation should be \geq 90%.

Optional: We routinely characterize our cultures at days 13, 16, and 25 to confirm the degree of vmDA specification and neuronal identity within the cultures. If desired, fix cultures or recover cells from wells at these stages for immunofluorescence and/or RT-qPCR. vmDA progenitors co-expressing LMX1A, EN1, NURR1 should be apparent by D16, and vmDA neurons co-expressing FOXA2 along with the postmitotic lineage markers PITX3 and TH should be readily apparent by D25 (Figure 2B; Table 5).

Α	В	С	D	Е	F	G	Н	I	J	K
	Nestin	OTX2	FOXA2	LMX1A	EN1	NURR1	PITX2	BARHL1	(PITX3) _e GFP	TH
	1:1000	1:1000	1:500	1:200	1:50	1:200	1:500	1:500	1:1000	1:1000
Day 13	+++	+++	+++	++	+/-	-/+	+/-	+/-	-	-
Day 16	+++	+++	+++	++	++	++	+/-	+/-	+/-	+/-
Day 25	++	+++	+++	++	+/-	+	+/-	+/-	++	++
Day 60	+/-	++	++	+	+/-	+/-	+/-	+/-	+++	+++

Table 5. Expression of vmDA Genes during Differentiation

From day 30 onwards, perform half-media changes, such that ~ **300 μL** is removed from each well and ~ **500 μL** of fresh, equilibrated media is added to maintain ~ **1 mL** /48-well. At these late stages of the differentiation, media changes can be done 2–3 times per week.

Note: For in vitro studies of DA development, homeostasis, disease modelling, or for drug screening approaches, vmDA neurons can be cultured until at least day 60, at which point electrophysiologically active vmDA neurons will be present.

Cryopreservation of vmDA Progenitors 23m

30 Timing: 1 h

This step describes the protocol for generating frozen stocks of vmDA progenitors for downstream applications.

Note: Cryopreserving vmDA progenitors between days 13 and 19 is ideal to balance lineage acquisition and survival.

- 31 Label cryovials and cool freezing container to $\, \, \& \, \, 4 \, \, ^{\circ} \text{C} \, \, .$
- Prepare freeze media containing [M]20 % DMSO and [M]60 % KnockOut serum replacement and [M]20 % fresh Maturation Media + BGCATD + [M]10 Micromolar (μM) Y27632.

Note: □0.5 mL freeze media and □0.5 mL Maturation Media are needed per 1 mL cryovial. Typically,

33 🔲 🧞 🎘

17m

Wash wells twice in PBS^{-/-} and add $\Box 400~\mu L$ /48-well of Accutase. Place in incubator for $\odot 00:07:00$ – $\odot 00:10:00$. Monitor cell morphology closely under the microscope.

34

Dissociate cells by gentle pipetting using a P1000 and dilute Accutase in Maturation Media + BGCATD + Y27632.

35

3m

Centrifuge cells at 300 x g, 00:03:00, aspirate supernatant and resuspend in Maturation Media + BGCATD + IM 10 Micromolar ($\textcircled{\mu M}$) Y27632 (1 mL/cm²) at approximately 1 mL x 10^6 /mL for accurate counting using a hemocytometer and Trypan blue exclusion.

36

3m

Centrifuge at $300 \times g$, 00:03:00, aspirate supernatant and resuspend in fresh Maturation Media + BGCATD + [M110 Micromolar (μ M) Y27632 at 6×10^6 /mL.

37

Transfer 300 μL of the cell solution to each cryovial and add 500 mL of freeze media, pipetting twice to gently mix.

39 **(II**)

2d

After $\geq \circlearrowleft \textbf{48:00:00}$, transfer cryovials to LN2 for long-term storage.

Pause Point: vmDA progenitors can be stored for periods ≥1 year.

CRITICAL: Ensure cells are exposed to Accutase for the minimum required time that still allows for accurate quantification. Maximizing the presence of small aggregates (2–10 cells) will improve cell survival. In contrast, over-dissociation in Accutase, leading to a high number of single cells, will significantly impair survival during thawing.

40 Timing: 30 min

See Figure 3 for an overview of vmDA progenitor maturation following thawing.

This step describes how to thaw cryopreserved vmDA progenitors for use in downstream applications.

41



Coat 48-well plates with [M]10 μ g/mL Laminin-521 and place in incubator for $\geq \bigcirc 02:00:00$.

42 Equilibrate \square 15 mL Maturation Media + BGCATD + [M]10 Micromolar (μ M) Y27632 at \emptyset 37 °C for \ge \bigcirc 00:15:00 .

43



Rapidly thaw a cryovial of vmDA progenitors in a § 37 °C water bath and transfer the contents to

■10 mL prewarmed media using a P1000 pipette.

44



3m

2h

Centrifuge at $300 \times g$, 00:03:00, aspirate supernatant and resuspend in Maturation Media + BGCATD + [M]10 Micromolar (μ M) Y2732.

45



Plate progenitors at $0.3-0.45 \times 10^6/\text{cm}^2$.

Optional: Count vmDA progenitors using a hemocytometer to ascertain an accurate count of surviving cells. Live/dead analysis using FACS can also be used with a subfraction of the cells, in addition to Trypan blue exclusion. Viability immediately post-thaw is typically 60%–70%, although a further proportion may undergo apoptosis after re-plating.

Optional: Transplant vmDA progenitors directly after thawing. Wash once in PBS^{-/-} after counting, centrifuge, aspirate supernatant and resuspend in Maturation Media + BGCATD + [μ]10 Micromolar (μ M) Y27632 at 0.1–0.2 x 10⁶/ μ L for transplantation.

Optional: For long-term culture of vmDA neurons for studying development, drug screening, and/or disease modeling, wash cells once in PBS^{-/-} the day after thawing and replace media with

■1 mL fresh Maturation Media + BGCATD. Proceed with steps 28–29 to continue the differentiation.

protocols.io

07/21/2021

Note: The concentration of rhGDNF and Y27632 can be doubled to [M]40 ng/mL and [M]20 Micromolar (μ M), respectively, if survival upon thawing is poor.

Immunostaining of vmDA Progenitors and Neurons

1h 42m

46 Timing: 2 days

This step outlines how to undertake immunocytochemical profiling of generated vmDA progenitors.

47

10m

Wash wells once with PBS- $^{\prime\prime}$ to remove debris and add $\ \Box 200 \ mL \ \ 4\% \ \ paraformal dehyde$. Leave at

§ Room temperature for © 00:10:00.

48 🔊 🕕

2m

Aspirate paraformaldehyde and wash 3 times with PBS^{-/-}, leaving wells in PBS^{-/-} for at least **© 00:02:00** between washes.

Pause Point: Fixed plates can be stored for several months at § 4 °C . Addition of [M]0.05 % Na- Azide to PBS^{-/-} can reduce possible bacterial or fungal contamination if storing long-term.

- 48.1 Aspirate paraformaldehyde.
- 48.2 Wash with PBS^{-/-}.
- 48.3 Leave wells in PBS $^{-/-}$ for at least \bigcirc 00:02:00 .

2m

48.4 Wash with PBS^{-/-}.

2m

48.5 Leave wells in PBS^{-/-} for at least \bigcirc **00:02:00**.

48	6	Wash with	DDC-/-
40.	()	wash with	PRS'

49 Prepare blocking buffer by adding 1 mL species-specific serum (same as secondary antibodies) and 0.6 mL 10% Triton X-100 to 18.4 mL PBS-Azide 0.05%.
 50 Add 200 μL /48-well of blocking buffer to each well for 00:30:00 - 01:00:00.

Dilute the primary antibodies in blocking buffer. Aspirate blocking buffer and add 200μ L /48- well of primary antibody solution. Incubate at $4 \, ^{\circ}$ C for 212:00:00 - 218:00:00.

52 A

Wash wells 3 times in PBS^{-/-} to remove primary antibodies and block cells again for \bigcirc **00:30:00** – \bigcirc **01:00:00** by adding \square **200** μ L blocking buffer .

- $52.1 \hspace{0.5cm} \text{Wash wells in PBS}^{\text{-/-}} \text{to remove primary antibodies}.$
- $52.2 \hspace{0.5cm} \text{Wash wells in PBS-}{}^{\prime\prime} \text{to remove primary antibodies}.$
- 52.3 Wash wells in PBS^{-/-} to remove primary antibodies.
- 52.4 Block cells again for \bigcirc 00:30:00 \bigcirc 01:00:00 by adding \square 200 μ L blocking buffer.
- 53 T

Dilute appropriate secondary antibodies in blocking buffer and add 200μ L secondary antibody solution to each well. Incubate cultures in the dark at 8000μ C Room temperature for 01:30:00 - 00:00:00.

54 Po

Wash wells in PBS^{-/-} and add [M]1 μ g/mL DAPI (diluted in PBS-/-) for \bigcirc 00:10:00 .



Wash wells twice in PBS^{-/-} and image.

Pause Point: Stained samples can be stored for a week at § 4 °C before imaging. It is best to store them in PBS^{-/-} + [M]**0.05 % Na-Azide**.

Note: It is optimal to perform immunofluorescence immediately after fixation of the cells. If long-term storage is necessary, ensure sufficient PBS-Azide is added to each well (🔄 1 mL), to account for evaporation during refrigeration and minimize exposure to light to prevent fluorochrome fading.