

NOV 11, 2022







MAP-2 & GFAP Staining Protocol

dx.doi.org/10.17504/protocols.io.q26g7y4b1gwz/v1

Machlusil Husna¹, Kusworini Handono², Hidayat Sujuti³, Aulanni'am⁴, Ettie Rukmigarsari⁵

- ¹Department of Neurology Faculty of Medicine, Universitas Brawijaya / dr. Saiful Anwar General Hospital, Malang, Indonesia;
- ²Department of Clinical Pathology Faculty of Medicine, Universitas Brawijaya / dr. Saiful Anwar General Hospital, Malang, Indonesia;
- ³Department of Ophthalmology Faculty of Medicine, Universitas Brawijaya / dr. Saiful Anwar General Hospital, Malang, Indonesia;
- ⁴Department of Chemistry, Faculty of Science Universitas Brawijaya, Malang, Indonesia;
- ⁵Faculty of Teacher Training and Education, Islamic University of Malang, Indonesia



ABSTRACT

Neurodegeneration due to neurotoxicity is one of the phenomena in temporal lobe epilepsy. Experimentally, hippocampal excitotoxicity process can occur due to kainic acid exposure, especially in the CA3 area. Neuronal death, astrocyte reactivity and increased calcium also occur in hippocampal excitotoxicity, but few studies have investigated immediate effect after kainic acid exposure. The organotypic hippocampal slice culture (OHSC) is a useful model for studying the neurodegeneration process, but there are still many protocol differences. In this study, minor modifications were made in the OHSC protocol

DO

dx.doi.org/10.17504/protocols.io.g26g7y4b1gwz/v1

PROTOCOL CITATION

Machlusil Husna, Kusworini Handono, Hidayat Sujuti, Aulanni'am, Ettie Rukmigarsari 2022. MAP-2 & GFAP Staining Protocol. **protocols.io**

https://dx.doi.org/10.17504/protocols.io.q26g7y4b1gwz/v1



-

Citation: Machlusil Husna, Kusworini Handono, Hidayat Sujuti, Aulanni'am, Ettie Rukmigarsari MAP-2 & GFAP Staining Protocol https://dx.doi.org/10.17504/protocols.io.g26g7y4b1gwz/v1

LICENSE

This is an open access protocol distributed under the terms of the <u>Creative Commons</u>

<u>Attribution License</u>, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

Oct 18, 2022

LAST MODIFIED

Nov 11, 2022

PROTOCOL INTEGER ID

71479

Sample Preparation

- 1 -Tissue slices were fixed using 4% PFA for 15 minutes
 - -Wash with PBS for 3 times (5 minutes each)
 - -Remove PBS and add 0,1% Triton X-100 for 30 minutes
 - -Wash with PBS for 3 times (5 minutes each)
 - -Add 10% BSA for 30 minutes at room temperature

Assay Procedure

- 2 -Incubate with primary antibody MAP-2 (1:500) and FDAP (1:500) overnight at temperature 4 °C.
 - -Wash with PBS for 3 times (5 minutes each)
 - -Incubate with secondary antibody for 1 h at room temperature
 - Goat anti-Rabbit IgG-F (1:500) for MAP-2 (Santa Cruz Biotechnology sc-53805, Lot#K1109)
 - Anti-mouse Rhodamine (1:500) for GFAP (Rockland, paint: 610-1002)
 - -Wash with PBS for 3 times (5 minutes each)
 - -Observe using CLSM with magnification of 100X and a wavelength of 488 nm (MAP-2) and 543 nm (GFAP)