



Sep 29, 2020

# Axy Prep DNA

1  
1

1 Works for me

This protocol is published without a DOI.

Zhang Xue



## PROTOCOL CITATION

Xue Zhang 2020. Axy Prep DNA. protocols.io  
<https://protocols.io/view/axy-prep-dna-bmuqk6vw>

## LICENSE

This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

## CREATED

Sep 29, 2020

## LAST MODIFIED

Sep 29, 2020

## PROTOCOL INTEGER ID

42608

1. Add 2ml EP 1200rcf 1min
2. Add 250µl Buffer S1 RNase A 4°C
3. Add 250µl Buffer S2 4-6 5min 00:05:00
4. Add 350µl Buffer S3 6-8 1200rcf 10min 00:10:00
5. Add 2ml 1200rcf 1min
6. Add 500µl Buffer W1 1200rcf 1min
7. Add 700µl Buffer W2 1200rcf 1min

8 离心

9 12000rcf 1min

10 1.5ml 40µl ddH<sub>2</sub>O 1-10min 12000rcf 1min ⌚ 00:05:00