



# SPARC\_Duke\_PelotGrill\_OT2-OD025340\_PigVagusNerve\_Morphology

Nicole A. Pelot<sup>1</sup>, Gabriel B. Goldhagen<sup>1</sup>, Jake E. Cariello<sup>1</sup>, Warren M. Grill<sup>1</sup>

<sup>1</sup>Duke University

1 Works for me [dx.doi.org/10.17504/protocols.io.6bvhan6](https://dx.doi.org/10.17504/protocols.io.6bvhan6)

## SPARC

Tech. support email: [info@neuinfo.org](mailto:info@neuinfo.org)



Nicole Pelot  
Duke University

## ABSTRACT

The protocol describes image segmentation and image analysis methods to quantify pig vagus nerve morphology from Masson's trichrome-stained nerve cross sections.

## DOI

[dx.doi.org/10.17504/protocols.io.6bvhan6](https://dx.doi.org/10.17504/protocols.io.6bvhan6)

## PROTOCOL CITATION

Nicole A. Pelot, Gabriel B. Goldhagen, Jake E. Cariello, Warren M. Grill 2020. SPARC\_Duke\_PelotGrill\_OT2-OD025340\_PigVagusNerve\_Morphology. **protocols.io**  
[dx.doi.org/10.17504/protocols.io.6bvhan6](https://dx.doi.org/10.17504/protocols.io.6bvhan6)

## KEYWORDS

Vagus nerve, nerve morphology, pig vagus nerve, endoneurium, perineurium, epineurium, fascicles, image segmentation

## LICENSE

————— This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

## CREATED

Aug 09, 2019

## LAST MODIFIED

Jun 30, 2020

## PROTOCOL INTEGER ID

26709

## MATERIALS TEXT

- Nikon's NIS Elements
- Matlab

## DISCLAIMER:

### DISCLAIMER – FOR INFORMATIONAL PURPOSES ONLY; USE AT YOUR OWN RISK

The protocol content here is for informational purposes only and does not constitute legal, medical, clinical, or safety advice, or otherwise; content added to [protocols.io](https://protocols.io) is not peer reviewed and may not have undergone a formal approval of any kind. Information presented in this protocol should not substitute for independent professional

## Image Segmentation

- 1 We used Nikon's NIS Elements software (v5.02.01, Build 1270) to segment pig vagus nerve micrographs stained with Masson's trichrome using the General Analysis RGB tool.
- 2 For each image, we selected preprocessing steps, such as smoothing and sharpening.
- 3 For each image, we selected ranges of hues, saturations, and intensities to values that identify the endoneurium (fascicles) and different values to identify the entirety of the nerve.
- 4 For each image, we selected postprocessing steps, such as setting a minimum size criterion (eliminate small off-target regions), smoothing, cleaning, closing, and filling holes.
- 5 We made manual adjustments as needed, including manual deletion of off-target regions and filling of target areas that had not been captured.
- 6 We converted the binary segmented image into "Graticule Masks", binary images saved as TIFs.

## Image Analysis

- 7 We imported the TIFs into Matlab and generated a data structure of the x and y coordinates of the pixels for each closed boundary of the loaded binary images using the *bwboundaries* function.
- 8 We scaled the pixel coordinates to microns using the segmented scale bar.
- 9 We calculated cross-sectional area of each fascicle and nerve using Matlab's polyarea. Effective diameter (for a nerve or fascicle) is the diameter of the circle that has the same cross-sectional area as the raw trace.