



DEC 04, 2022

WORKS FOR ME

1

HuBMAP Lung TMC UCSD & URM Protocol Overall for Multimodal Single Cell/Nucleus Assays

COMMENTS 0

DOI

dx.doi.org/10.17504/protocols.io.4r3l27qpqg1y/v1

[Gloria S Pryhuber](#)¹, [Elizabeth Duong](#)²,
[Blue Lake](#)², [Cory Poole](#)¹, [heidie_huyck](#)¹,
[Ravi Misra](#)¹, [Xin Sun](#)², [Jim Hagood](#)³,
[Dinh H Diep](#)², [Kun Zhang](#)²

¹University of Rochester Medical Center;²University of California San Diego;³University of North Carolina at Chapel Hill

Human BioMolecular Atlas Program (HuBMAP)
Method Development Community

LungMap2 Consortium

[1 mc](#)

Gloria S Pryhuber

[University of Rochester Medical Center](#)

ABSTRACT

This protocol is an overview of the protocols in use for the single nuclear RNA sequencing and ATACseq datasets for HubMAP-Lung Tissue Mapping Center, TMC UCSD / TMC URM. It includes protocols.io links to each of the individual protocols that make up this project workflow, from processing tissue, sectioning, H&E imaging the tissue samples, isolation of nuclei from frozen lung tissue, running of multimodal single cell/nucleus assays.uploading the final data to the HuBMAP HIVE.

See also: Human Kidney, Urinary Tract, and Lung Cell Type Mapping Pipeline for the Human Biomolecular Atlas Program (HuBMAP) V.2 dx.doi.org/10.17504/protocols.io.bvvin64e

DOI

dx.doi.org/10.17504/protocols.io.4r3l27qpqg1y/v1

PROTOCOL CITATION

Gloria S Pryhuber, Elizabeth Duong, Blue Lake, Cory Poole, heidie_huyck, Ravi Misra, Xin Sun, Jim Hagood, Dinh H Diep, Kun Zhang 2022. HuBMAP Lung TMC UCSD & URM Protocol Overall for Multimodal Single Cell/Nucleus Assays . **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.4r3l27qpqg1y/v1>



FUNDERS ACKNOWLEDGEMENT

HuBMAP

Grant ID: U54 HL145608; U54HL165443

LICENSE

This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

Dec 04, 2022

LAST MODIFIED

Dec 04, 2022

PROTOCOL INTEGER ID

73519

MATERIALS TEXT

Required materials are provided within each linked protocol.

- 1 Screen and accept donor organs, gross and preserve the tissues into samples for analysis using the protocols:

602.2 Donor Acceptance Criteria for URMHC HTC HuBMAP and LungMAP Inclusion
dx.doi.org/10.17504/protocols.io.bjuxknxn
603.3 & 604.5_URMC_HTC_Whole Lung and Lobe Processing
dx.doi.org/10.17504/protocols.io.biz7kf9n
611.2 URMHC HTC Formalin-Inflated, Paraffin-Embedded Human Lung Tissue
dx.doi.org/10.17504/protocols.io.bjttknkn
615.1 URMHC HTC Non-Inflated Fresh-Frozen Embedded Lung and Associated Tissue
dx.doi.org/10.17504/protocols.io.bjtnknme
615.2 URMHC HTC Non-Inflated Fresh-Frozen Embedded Lung and Associated Tissue V2
dx.doi.org/10.17504/protocols.io.rm7vz84k2vx1/v2
- 2 Create the cryosections to be used in the single nuclear isolation preparation. Generally placing 10 x 40 micron sections in a single cryovial.
624.1_URMC_HTC_Cryosectioning Non-Inflated Fresh-Frozen Embedded Lung for Nuclear Isolation or Multiplexed In Situ Hybridization
dx.doi.org/10.17504/protocols.io.3byl4j91jlo5/v1
- 3 Register via the HuBMAP UUID generator at <https://ingest.hubmapconsortium.org/>, the Donor, the Organ Lobe, the Tissue Block, and the set of tissue sections (40 micron curls) from that block used in single nuclear isolation. common coordinate region information. Record and keep locally the key for the generated UUIDs.
- 4 Register the Tissue Blocks, just registered in the UUID generator, in the Common Coordinate Framework via the Registration User Interface at <https://hubmapconsortium.github.io/ccf-ui/rui/>
- 5 Approximately every 10th section, make a 10 micron section just before and just after, the set of 40 micron curls. Stain these thinner sections by H&E.
621.1.HTC_H&E Stain (Paraffin or Cryosections) dx.doi.org/10.17504/protocols.io.36wggjinq3vk5/v1
- 6 Isolate single nuclei from the tissue curls using protocol:
Frozen lung tissue dissociation for single-nucleus RNA & Chromatin assays
dx.doi.org/10.17504/protocols.io.bh26j8he

- 7 Running of multimodal single cell/nucleus assays.
10X snRNA-seq: dx.doi.org/10.17504/protocols.io.86khzcw
10X snATAC-seq: dx.doi.org/10.17504/protocols.io.bvssn6ee
SNARE-Seq2: dx.doi.org/10.17504/protocols.io.bzdsp26e

- 8 To Be Added:
10X Chromium Single Cell *Multiome* ATAC + Gene Expression
MERFISH
DARTFISH