




APR 15, 2024

## In vivo Ephys: Analysis Protocol

 In 2 collections

Sasha Burwell<sup>1</sup>

<sup>1</sup>Duke University

ASAP Collaborative Research Network



Sasha Burwell  
Duke University

OPEN  ACCESS



### ABSTRACT

This protocol details the analysis procedures performed on the in vivo electrophysiology data post-spike sorting.

#### DOI:

[dx.doi.org/10.17504/protocols.io.261ge5o2yg47/v1](https://dx.doi.org/10.17504/protocols.io.261ge5o2yg47/v1)

**Protocol Citation:** Sasha Burwell 2024. In vivo Ephys: Analysis Protocol. **protocols.io**  
<https://dx.doi.org/10.17504/protocols.io.261ge5o2yg47/v1>

**License:** This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working

**Created:** Apr 02, 2024

**Last Modified:** Apr 15, 2024

## Analysis Protocol

- 1 To plot the waveforms from one electrode of interest in the raw data (ex, Fig. 1c), run the Matlab code

```
plotRawTraces(electrode)
```

on the raw continuous data.

- 2 To analyze the sorted templates extracted from Spyking Circus,

- 2.1 Use the custom matlab code.

```
plotSpikes(chWithCell)
```

- 2.2 chWithCell = array of the electrode which has a cell for each template (ex [5 15 3]).

### Note

Make sure this is in the order of the templates output by Spyking Circus.

- 2.3 Calculates and plots spike amplitude, cumulative spikes, and firing rate in Hz for all the templates saved.

- 3 From the resulting Templates\_data structure, save each row as an individual dat variable.

i.e., for ch 15 in the above example:

```
load Templates_data;
dat = Templates_data(2, :); %row 2 is channel 15 data
save('mouse_date_ch15_data', 'dat');
```

#### 4 Copy multiple channels into one folder for grouped analysis, and run

```
spikingAcrossCells
```

in this folder.

#### 5 Referencing the control and experimental

```
allCells
```

data structures created from

```
spikingAcrossCells
```

, run

```
plotFeaturesAcrossTime
```

to generate the plots seen in the left of Fig. 1d, e, and f.

#### 6 Use the features in the

```
allCells
```

data structure:

**6.1** Calculate delta norm as  $(\text{post-pre})/(\text{post+pre})$  for each feature for each cell.

**6.2** Input ddHTP values as

```
y1
```

and +HTP values as

```
y2
```

into the Matlab code

```
TwoGroupsMeanDiff
```

for statistical analysis and generating the plots seen in the right of Fig. 1d, e, and f.

- 6.3** Use GraphPad Prism to generate the plots in Ext. Fig 2d and h-l from the delta norm features.