



FEB 13, 2024

🌐 Catalepsy test in rats

eduard.bentea¹, María Sanchiz Calvo¹, Veerle Baekelandt¹

¹Katholieke Universiteit Leuven



eduard.bentea

ABSTRACT

Protocol for performing the catalepsy bar test in rats. This test evaluates muscular rigidity, and can be used in rat models with motor deficits, such as models of Parkinson's disease.

OPEN  ACCESS



DOI:

dx.doi.org/10.17504/protocols.io.q26g7p24kgwz/v1

Protocol Citation: eduard.bentea, María Sanchiz Calvo, Veerle Baekelandt 2024. Catalepsy test in rats. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.q26g7p24kgwz/v1>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working
We use this protocol and it's working

Created: Feb 13, 2024

Last Modified: Feb 13, 2024


PROTOCOL integer ID: 95155

Keywords: ASAPCRN

Funders Acknowledgement:
ASAP (Aligning Science Across
Parkinson's)

Test

1h

- 1 Bring cages to the behavioral room for at least  01:00:00 prior to the test for habituation

1h
- 2 Place both forepaws of the rat on an horizontal bar elevated so that the rat is in an half rearing position
- 3 Measure the descent latency, defined as the time it takes for both forepaws to be removed from the bar
- 4 Repeat the test three times, with a brief rest in between