



Nov 02, 2021

Virus Production and Administration

Sarkis Mazmanian¹, Bryan Yoo¹, Jessica Griffiths¹

¹California Institute of Technology

1

dx.doi.org/10.17504/protocols.io.bzn6p5he

Mazmanian Lab

Jessica Griffiths

Protocol for virus production and administration to mice in Yoo et al 2021

DOI

dx.doi.org/10.17504/protocols.io.bzn6p5he

Sarkis Mazmanian, Bryan Yoo, Jessica Griffiths 2021. Virus Production and Administration. **protocols.io**

https://dx.doi.org/10.17504/protocols.io.bzn6p5he

Emerald Foundation

Grant ID: n/a

National Institutes of Health

Grant ID: GM007616 and DK078938

Department of Defense

Grant ID: PD160030

Center for Microbial and Environmental Interactions (CEMI)

Grant ID: n/a

Aligning Science Across Parkinson's

Grant ID: ASAP-000375

 $_{\scriptscriptstyle
m oldsymbol{ iny}}$ protocol ,

Nov 02, 2021

Nov 02, 2021



1

Virus Production

- 1 Virus was produced using methods described in Challis et. al. (Challis et al., 2019)
- 2 Cells were grown in DMEM + Glutamax + Pyruvate (Gibco, Gaithersburg, MD-Stock# 10569-010) + 5% FBS + non-essential amino acids (Gibco, Gaithersburg, MD-Stock# 11140-050), penicillin-streptomycin (Gibco, Gaithersburg, MD-Stock# 15070-063).
- 3 Briefly, human embryonic kidney (HEK293T) cells were triple-transfected with pUCmini-iCAP-AAV-PHP.S, pHelper Plasmid, and one of the following pAAV genomes (CAG-NLS-GFP, hSYN1-tDTomato, GFAP-tDTomato, hSYN1-mRuby2, hSYN1-DiO-mRuby2, hSYN1-mNeonGreen, hSYN1-DiO-mNeonGreen, hSYN1-mTurquoise2, hSYN1-DiO-mTurquoise2, hSYN1-DiO-hM3Dq-mRuby2, CAG-GCaMP6F) in a DPBS + polyethylenimine (PEI)
- 4 Virus is precipitated from cells and supernatant with an 8% PEG solution (wt/vol), and purified by ultracentrifugation using 15%, 25%, 40%, 60% stacked iodixanol gradients.

Virus Administration

- 5 Virus was titered to 10^{12} viral genomes (vgs) and resuspended to a volume of 100μ l with sterile PBS.
- 6 Mice were anesthetized using 2% isoflurane, and virus was injected retro-orbitally.