



OCT 18, 2023

OPEN ACCESS



**DOI:**  
[dx.doi.org/10.17504/protocols.io.81wgbx55qlpk/v1](https://dx.doi.org/10.17504/protocols.io.81wgbx55qlpk/v1)

**Protocol Citation:** Núria Peñuelas 2023. Blood pressure. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.81wgbx55qlpk/v1>

**License:** This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working  
We use this protocol and it's working

**Created:** Oct 18, 2023

**Last Modified:** Oct 18, 2023


## Blood pressure

Núria  
Peñuelas<sup>1</sup>

<sup>1</sup>Vall d'Hebron Research Institute

Vilalab Public

Vilalab Public



Núria Peñuelas  
Vall d'Hebron Research Institute

### ABSTRACT

Blood pressure measurement in mice

- 1** Measure blood pressure using the tail-cuff method (LE-5002 Non-Invasive Blood Pressure Meter, Panlab).
- 2** Assess the systolic and diastolic blood pressure (in mmHg) in two independent sessions for each mouse.
- 3** Calculate the mean from the valid values of 10 measurements in each session.