


JUN 01, 2023


Cremophor EL agar (EL slant)

 In 1 collection

Amy

Gladfelter¹

¹Duke University



Amy Gladfelter

ABSTRACT

This media is used to determine the ability of *Malassezia* species to utilize Cremophor EL.

ATTACHMENTS

[Malassezia media_1.pdf](#)

GUIDELINES

This media is used to determine the ability of *Malassezia* species to utilize Cremophor EL.

OPEN ACCESS

DOI:

dx.doi.org/10.17504/protocols.io.n92ld3m87g5b/v1

Protocol Citation: Amy Gladfelter 2023. Cremophor EL agar (EL slant).

protocols.io

<https://dx.doi.org/10.17504/protocols.io.n92ld3m87g5b/v1>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working
We use this protocol and it's working


Created: Nov 03, 2018

Last Modified: Jun 01, 2023

PROTOCOL integer ID:

17459

1 SDA  65 g

2 Cremophor EL (Sigma)  10 mL

3 dH₂O  900 mL

4 Sterilize by autoclaving.