

Apr 18, 2024

# 10x Protocols: Chromium Next GEM Single Cell 3' -- University of Minnesota TMCs (CG000315 Rev E)

DOI

[dx.doi.org/10.17504/protocols.io.261ge5ozdg47/v1](https://dx.doi.org/10.17504/protocols.io.261ge5ozdg47/v1)

Laura J Niedernhofer<sup>1</sup>, David A Bernlohr<sup>1</sup>

<sup>1</sup>University of Minnesota, Minneapolis, MN

Cellular Senescence Net...

UMN SenNet



Allie

Allie Pybas

UMN

OPEN  ACCESS



DOI: [dx.doi.org/10.17504/protocols.io.261ge5ozdg47/v1](https://dx.doi.org/10.17504/protocols.io.261ge5ozdg47/v1)

**Protocol Citation:** Laura J Niedernhofer, David A Bernlohr 2024. 10x Protocols: Chromium Next GEM Single Cell 3' -- University of Minnesota TMCs (CG000315 Rev E). protocols.io <https://dx.doi.org/10.17504/protocols.io.261ge5ozdg47/v1>

**License:** This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working

**We use this protocol and it's working**

**Created:** April 18, 2024

**Last Modified:** April 18, 2024

**Protocol Integer ID:** 98430

**Funders Acknowledgement:**

NIH

Grant ID: 5U54AG076041-03

NIH

Grant ID: 5U54AG079754-02



## Abstract

Following single-cell dissociation or nuclei isolation, the Next GEM 3' assay uses microfluidics to partition and assign cell- or nuclei-specific barcodes to transcript cDNA at the 3' end. Barcoded cDNA is prepared with adaptors for sequencing by synthesis (SBS). The following protocol has been adapted from protocols developed by 10x Genomics and Illumina to be used at the University of Minnesota TMCs in collaboration with the University of Minnesota Genomics Center. These protocols are owned by their respective companies and are subject to periodic revision.



## Tissue Preparation

- 1 Complete single cell or nuclei isolation prior to starting this protocol

## Library Preparation

- 2  CG000315\_Rev-E (4).pdf 4.2MB

### Note

Sequence with the read format 29,10,10,89

### Note

Sequencers used at UMN Genomics Center:

- Illumina NextSeq 2000
- Illumina NovaSeq 6000
- Illumina NovaSeq X Plus

## FASTQ Generation

- 3 BCL data from Illumina sequencer is demultiplexed and converted into FASTQ format using bcl2fastq version 2.20.0