

DEC 08, 2023

## UroMOCA and ColoMOCA implantation protocol


Brett

Hanzlicek<sup>1</sup>, Margot Damsaser<sup>2</sup>, Dennis Bourbeau<sup>3</sup>

<sup>1</sup>Advanced Platform Technology Center, Louis Stokes VA Hospital, Cleveland, Ohio;

<sup>2</sup>Department of Biomedical Engineering, Cleveland Clinic Lerner Research Institute, Cleveland, Ohio;

<sup>3</sup>MetroHealth Medical Center, Cleveland, Ohio



Brett Hanzlicek

### ABSTRACT

This protocol describes the implantation of our bladder (UroMOCA) and bowel device (ColoMOCA) in pigs.

OPEN  ACCESS



### DOI:

[dx.doi.org/10.17504/protocols.io.rm7vzxq78gx1/v1](https://dx.doi.org/10.17504/protocols.io.rm7vzxq78gx1/v1)

**Protocol Citation:** Brett Hanzlicek, Margot Damsaser, Dennis Bourbeau 2023. UroMOCA and ColoMOCA implantation protocol .

**protocols.io**

<https://dx.doi.org/10.17504/protocols.io.rm7vzxq78gx1/v1>

**License:** This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working  
We use this protocol and it's working

**Created:** Dec 05, 2023

## Anesthesia

- 1 Anesthetize/prep the animal for surgery.

Aseptic procedures by trained personnel are utilized during the surgery. All tools are sterilized prior to insertion. Drapes will be used and personnel performing surgery will wear clean surgical scrubs, sterile gowns/gloves and masks.

- 2 Insert single lumen 8 Fr catheter, and drain bladder. Save some urine for urinalysis.

## Laparotomy

- 3 Place animal in supine position.

- 4 Administer analgesic at injection site.

- 5 Make a midline suprapubic incision ~9 cm long to expose the bladder, peritoneal cavity and viscera.

- 6 Isolate bladder and incise the bladder dome.

## UroMOCA placement and bladder closure

- 7 Place UroMOCA device into bladder.

8 Close bladder incision and evaluate integrity.

9 Place sterile drape on abdomen for UroMOCA transmission coil.

10 Verify function of UroMOCA in vivo.

## Colotomy

11 Isolate the region of the colon where the implant will go.

12 Make a single, small incision into the colon.

## Insert and anchor ColoMOCA

13 Insert device through the small incision and secure in place with sterile non-absorbable sutures.

## Test Wireless Transmission

14 Wake the coloMOCA with radio if not already awake. Confirm data transmission.

## Close incision

- 15
- Close the colon with an absorbable suture in an interrupted pattern
  - Close the peritoneum and fascia with absorbable sutures using a continuous suture pattern.
  - Close the subcutaneous tissue and skin with non-absorbable monofilament suture with a cutting needle, using an interrupted suture pattern. Skin staples may also be used.

## Test Wireless Transmission (2)

- 16
- Confirm data transmission from both the UroMOCA and ColoMOCA, and then put devices to sleep

## CT-Imaging

- 17
- Image to verify placement of the sensors.

## Place jacket on the pig

- 18
- Place and zip the jacket on the pig while still anesthetized.

## Wake up the animal

- 19
- Stay with the animal until fully awake. Return to pig room.