


Apr 17, 2024

## Quality control assessment for microbial genomes: GalaxyTrakr MicroRunQC workflow

 Forked from [Quality control assessment for microbial genomes: GalaxyTrakr MicroRunQC workflow](#)

DOI

[dx.doi.org/10.17504/protocols.io.4r3l27oz3g1y/v1](https://dx.doi.org/10.17504/protocols.io.4r3l27oz3g1y/v1)



Ruth Timme<sup>1</sup>, Sai Laxmi Gubbala Venkata<sup>2</sup>, Maria Balkey<sup>1</sup>, Robyn Randolph<sup>3</sup>, William Wolfgang<sup>4</sup>, Errol Strain<sup>3</sup>

<sup>1</sup>US Food and Drug Administration; <sup>2</sup>Bacteriology Laboratory, Wadsworth Center, New York State Department of Health, Albany, New York, USA;

<sup>3</sup>Center for Veterinary Medicine, U.S. Food and Drug Administration, College Park, Maryland, USA; <sup>4</sup>Wadsworth Center NYSDOH

GenomeTrakr  
 Tech. support email: [genomeTrakr@fda.hhs.gov](mailto:genomeTrakr@fda.hhs.gov)

 **Maria Balkey**  
 US Food and Drug Administration

OPEN  ACCESS



DOI: [dx.doi.org/10.17504/protocols.io.4r3l27oz3g1y/v1](https://dx.doi.org/10.17504/protocols.io.4r3l27oz3g1y/v1)

External link: <https://galaxytrakr.org>

**Protocol Citation:** Ruth Timme, Sai Laxmi Gubbala Venkata, Maria Balkey, Robyn Randolph, William Wolfgang, Errol Strain 2024. Quality control assessment for microbial genomes: GalaxyTrakr MicroRunQC workflow. [protocols.io](#) <https://dx.doi.org/10.17504/protocols.io.4r3l27oz3g1y/v1>

**Manuscript citation:**

Timme, R. E., W. J. Wolfgang, M. Balkey, S. L. G. Venkata, R. Randolph, M. Allard, and E. Strain. 2020. Optimizing open data to support one health: best practices to ensure interoperability of genomic data from bacterial pathogens. *One Health Outlook* 2: 20.

**License:** This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working

**Created:** October 26, 2022

**Last Modified:** April 17, 2024

**Protocol Integer ID:** 71843

**Keywords:** WGS, Quality Control, GalaxyTrakr, GenomeTrakr, microbial pathogen surveillance

### Disclaimer

Please note that this protocol is public domain, which supersedes the CC-BY license default used by protocols.io.



## Abstract

**PURPOSE:** Step-by-step instructions for checking WGS sequence quality. The MicroRunQC workflow, implemented in a custom Galaxy instance, will produce quality assessments for raw reads (Illumina paired-end fastq files) and draft de novo assemblies, along with reporting the sequence type for each isolate. This workflow will work on most microbial pathogens, so we advise laboratories to upload their entire MiSeq/NextSeq run through this workflow.

**SCOPE:** This protocol covers the following tasks:

1. set up an account in GalaxyTrakr
2. Create a new history/workspace
3. Upload data
4. Execute the MicroRunQC workflow
5. Interpret the results

V3: updated with *Cronobacter* thresholds

## Account set up

- 1 1. Create a GalaxyTrakr account here: <https://account.galaxytrakr.org/Account/Register>

### User Registration Form

Location

Add New Location

First Name

Last Name

Email

Email will be used for automated messages to include registration information!

Primary Phone

If possible please use a mobile number than can accept text messages, only used for support!

Title

Requirements

- 1.1 Log into your GalaxyTrakr account: <https://galaxytrakr.org>

**Galaxy / GalaxyTrakr 1905** Analyze Data Workflow Visualize Shared Data Help Login

Welcome to Galaxy, please log in

Username or Email Address

Password

Forgot password? Click here to reset your password.

Don't have an account? Registration for this Galaxy instance is disabled. Please contact an administrator for assistance.

**Welcome to GalaxyTrakr: open-source bioinformatics for public health.**

**This site is intended for use by GenomeTrakr laboratories and their collaborators to assist in the analysis of genomic data for foodborne pathogens. This instance of Galaxy is hosted in a public environment and no personally identifiable (PII) or commercial confidential information should be uploaded.**

**---!!--Information and Announcements---!!--**

**Please re-import the skseamist workflow that was updated a few days ago. Previous versions are no longer working and are causing errors when running. Thank you.**

Access CFSAN SNP Pipeline workflows in the shared workflows screen.

Post in the official Galaxy GenomeTrakr board on the Redmine Site: [Click here](#)

[Click here](#) to access the GalaxyTrakr User Guide

Forgot Password? Email [GalaxyTrakr@SupportTeam](mailto:GalaxyTrakr@SupportTeam)

## Create a new history

- 2 **Create a new history.**

We recommend creating a new history for each new MiSeq Run and including the flow-cell ID and date in the history name.

Save your MicroRunQC output here and any other relevant analyses, like serotyping, or AMR detection.

After all the analysis output from this run is saved to your internal data network or computer, older history's should be purged/deleted so as not to occupy the limited storage space in your account. In some cases it may be useful to save, for a limited time, multiple histories or to run analyses concurrently in multiple histories. In these cases you need to pay attention to your % usage bar (shows % used of allocated storage space) in the upper right corner of the GalaxyTrakr page. If you need additional space you can contact [galaxytrakrsupport@fda.hhs.gov](mailto:galaxytrakrsupport@fda.hhs.gov) and request additional storage.

- 2.1 Click on the + icon in the upper right History panel

**History**

**test-runQC**

18 shown, 12 deleted, 22 hidden

750 b ☐ ☐ ☐




- 2.2 Name your new History by clicking on the "Unnamed history", type in desired name and hit enter. We recommend including the run cell ID and the date the run was started.




## Upload data

- 3 **This section will describe the process for uploading raw fastq files into your active History panel.** After the files have been uploaded they will stay in your account until they are deleted.

- 3.1 Click on the upload/download icon on the top of the left web page to start an upload process.



- 3.2 Select "Type (set all):auto-detect." Choose local file button and navigate to the desired fastq files, then click "start" to upload files. These files should be paired (two per sample/isolate).



As the file uploads complete, each row will turn green. Samples in yellow are still in process.

- 3.3 You have just upload a set of forward and reverse reads. For further analysis these files need to be paired properly so the platform knows which R1 and R2 files go with each sample/isolate. GalaxyTrakr does this by creating a **List of Dataset Pairs**.

Within your newly created History panel, click the "check box," then select all the files you just uploaded by clicking "All" or by individually selecting the ones you want to pair.




**Screenshot of History panel showing recently uploaded files.** Note the way the files are named, using R1 and R2 to identify the paired reads. This will be important in the next step. Some naming conventions can be slightly different.

3.4 Click "For all selected" and choose "Build List of Dataset Pairs"



3.5 A new window will open to help you pair the fastq files properly. Note how your paired reads are named (\_R1 and \_R2 in the example above)

Select **Clear filters**, then click **Auto-pair**.



Alternatively, instead of autopairing you can click "choose filters" and select the appropriate filter for the pairing:

### 3.6 Paired reads will pair in the middle column and turn green.




Name your dataset: Example, "pairedSet-<FlowCell>-<date>"


Click **Create list**.

Create a collection of paired datasets

### 3.7 This paired dataset will now be available for analysis in your history panel. You can run multiple analyses on the same dataset in a history rather than upload the same sequence data to a new history to perform additional analyses. This will help you use your allocated storage space efficiently.





**History**   

search datasets 

**test - SOP dev**

5 shown, 2 deleted

658.07 MB  


All None For all selected...


- ☐ 7: PairedList-RunCell-Date  
a list of pairs with 2 items
- ☐ 4: CFSAN074384\_S20\_L001\_R2\_00  
1.fastq.gz
- ☐ 3: CFSAN074384\_S20\_L001\_R1\_00  
1.fastq.gz
- ☐ 2: CFSAN074382\_S15\_L001\_R2\_00  
1.fastq.gz
- ☐ 1: CFSAN074382\_S15\_L001\_R1\_001.  
fastq.gz


Run the MicroRunQC workflow

4 **Add the MicroRunQC workflow to your own "workflows" panel.** You only have to do this step once for each new workflow you need.

4.1 Navigate to the "Shared Data" drop down menu, choose workflows and from the MicroRunQC drop down menu select import.

**Galaxy / GalaxyTrakr 1901** Analyze Data Workflow Visualize Shared Data Help User 

**Tools** 

search tools 

**Get Data**

NGS TOOLBOX


- NGS: QC and manipulation
- NGS: Assembly
- NGS: Screening and Prediction
- NGS: Mapping
- NGS: Annotations
- NGS: CFSAN SNP Pipeline (beta)
- NGS: Phylogenetics

METAGENOMICS TOOLBOX

- Metagenomics: Metaphlan
- Metagenomics: Diamond
- Metagenomics: Kraken
- Metagenomics: Mitokmer
- Metagenomics: Lefse
- Metagenomics: Assembly

UTILITIES


**Published Workflows**

search name, annotation, owner 

Advanced Search

Name	Annotation	Owner	Pages	Immunity Tags	Last Updated
MicroRunKraken2		estrain	★★★★★		Jul 23, 2019
MicroRunQC (formerly skesaMLST)		estrain	★★★★★		Jul 23, 2019
srr22_paired_list	SRR22 workflow for paired lists	estrain	★★★★★		May 28, 2019
CFSAN SNP Pipeline w/ Assembly Lite		jpayne	★★★★★		Mar 26, 2019
QC_Reads	FASTQC Run on List of Paired MiSeq FASTQ Files	estrain	★★★★★		Nov 02, 2018
CFSAN SNP Pipeline w/ SRA Download and Assembly		jpayne	★★★★★		Oct 30, 2018
CFSAN SNP Pipeline		jpayne	★★★★★		Oct 30, 2018

**Data Libraries:**

- Histories
- Workflows 
- Visualizations
- Pages
- Immunity Tags
- Last Updated

4.2 To see the new workflow in your "Workflows" tools panel on the left, open the Workflow tab and check "show in tools panel" for the workflow of interest.



Galaxy / GalaxyTrakr 1901

Analyze Data **Workflow** Visualize Shared Data Help User

Tools

search tools

Get Data

NGS TOOLBOX

NGS: QC and manipulation

NGS: Assembly

NGS: Screening and Prediction

NGS: Mapping

NGS: Annotations

NGS: CFSAN SNP Pipeline (beta)

NGS: Phylogenetics

METAGENOMICS TOOLBOX

Metagenomics: Metaphlan

Metagenomics: Diamond

Metagenomics: Kraken

Metagenomics: Mitokmer

Metagenomics: Lefse

Metagenomics: Assembly

UTILITIES

Your workflows

search for workflow...

Name	Tags	Owner	# of Steps	Published	Show in tools panel
Imported: MicroRunKraken2		You	8	No	<input checked="" type="checkbox"/>
Imported: MicroRunQC (formerly skesaMLST)		You	4	No	<input checked="" type="checkbox"/>
Imported: QC_Reads		You	5	No	<input checked="" type="checkbox"/>
Imported: CFSAN SNP Pipeline w/ SRA Download and Assembly		You	16	No	<input checked="" type="checkbox"/>
Imported: CFSAN SNP Pipeline w/ Assembly		You	16	No	<input checked="" type="checkbox"/>
Imported: skesaMLST		You	4	No	<input type="checkbox"/>
Imported: SeqSero by List of SRR Numbers		You	3	No	<input checked="" type="checkbox"/>
Imported: CFSAN SNP Pipeline (Lite)		You	9	No	<input type="checkbox"/>
Imported: CFSAN SNP Pipeline		You	13	No	<input checked="" type="checkbox"/>

#### 4.3 From the workflow menus select MicroRunQC

Galaxy / GalaxyTrakr 1901

Tools

search tools

NGS: Screening and Prediction

NGS: Mapping

NGS: Annotations

NGS: CFSAN SNP Pipeline (beta)

NGS: Phylogenetics

METAGENOMICS TOOLBOX

Metagenomics: Metaphlan

Metagenomics: Diamond

Metagenomics: Kraken

Metagenomics: Mitokmer

Metagenomics: Lefse

Metagenomics: Assembly

UTILITIES

Text Manipulation

Join, Subtract and Group

Graph/Display Data

Convert Formats

Filter and Sort

Workflows

All workflows

Imported: CFSAN SNP Pipeline

Imported: SeqSero by List of SRR Numbers

Imported: CFSAN SNP Pipeline w/ Assembly

Imported: CFSAN SNP Pipeline w/ SRA Download and Assembly

Imported: MicroRunKraken2


Imported: MicroRunQC (formerly skesaMLST)

Imported: QC\_Reads

#### 4.4 Select paired list dataset you created earlier.


Click Run Workflow. This can take some time depending on the number of samples you are analyzing. If you choose to you can log out of GalaxyTrakr and log back in at a later time to see if the job is completed.





4.5 Upon completion of the pipeline all tiles in the history bar will be green.


In the "Filter on Data" tile click on the "Eye" icon to view the output in the GalaxyTrakr window.



## Interpret the results

### 5 Download and interpret the results:

5.1 Click "Filter on data" and then the floppy disc icon. The tabular file can be opened in a text reader or converted to a format that can be opened on excel.



5.2 The MicroRunQC output file includes the following metrics:

Parameter	Input	Description
Contigs	Assembly	Number of contigs in the de-novo SKESA assembly. Contigs smaller than 200 base-pairs (bp) are not counted.



<b>Length</b>	Assembly	Total length of all contigs > 200bp. This should approximate the size of the genome for the target organism.
<b>EstCov</b>	Assembly	Mean coverage for contigs in the SKESA assembly.
<b>N50</b>	Assembly	Sequence length of the shortest contig at 50% of the total genome length
<b>MedianInsert</b>	Read	Distance between forward and reverse reads. Calculated by mapping reads to SKESA assembly using bwa.
<b>MeanLength_R1</b>	Read	Mean length of forward read
<b>MeanLength_R2</b>	Read	Mean length of reverse read
<b>MeanQ_R1</b>	Read	Mean Q-score of forward read
<b>MeanQ_R2</b>	Read	Mean Q-score of reverse read
<b>Scheme</b>	Assembly	PubMLST (pubmlst.org) database scheme (e.g. enterica for <i>Salmonella enterica</i> )
<b>ST</b>	Assembly	Sequence Type
<b>Loci</b>	Assembly	gene (allele number) – for example aroC(118)

**MicroRunQC output table headers.** This table lists the summary metrics for sequence quality, number of contigs, and estimated genome size, along with other common metrics for reads (Median Insert Size and Mean Length) and assemblies (N50). Additionally, if the Multi-Locus Sequence Type (MLST) for the isolate is available from pubmlst, the workflow also reports Sequence Type (ST) and the associated alleles.

\*\*This output should be saved either to your LIMS or to a spreadsheet linked to the sequencing run and samples.

### 5.3 Example output for 4 *Listeria* samples run through the MicroRunQC workflow:

File name	Contigs	Length	EstCov	N50	MedianInsert	MeanLength_R1	MeanLength_R2	MeanQ_R1	MeanQ_R2	Scheme	ST				
FSL-R9-8346	14	2876874	87.8	512255	408	147.6	147.7	33.1	32.7	lmonocytogenes	389	abcZ(52)	bglA(1)	cat(12)	dap
FSL-R9-8348	11	2832172	84.2	1464158	388	147.9	147.9	33.2	32.9	lmonocytogenes	795	abcZ(7)	bglA(10)	cat(18)	dap
FSL-R9-8350	14	2884629	64.2	450082	390	147.2	147.2	33.1	32.7	lmonocytogenes	37	abcZ(5)	bglA(7)	cat(3)	dap
FSL-R9-8352	12	2902520	85.9	1460419	390	148.1	148.2	33.1	32.8	lmonocytogenes	391	abcZ(7)	bglA(6)	cat(62)	dap

Spreadsheet showing example output for 5 *Listeria monocytogenes* samples from a NextSeq sequencing run.

### 5.4 Quality control threshold guidelines for the GenomeTrakr surveillance network. These are also relevant for NARMS and VetLIRN contributors.

\*MicroRunQC users should follow threshold guidelines established by their respective surveillance coordinating body(s).

A	B	C	D	E	F	G	H
Quality metric	<b><i>Salmonella</i></b>	<b><i>Listeria</i></b>	<b><i>E. coli</i></b>	<b><i>Shigella</i></b>	<b><i>Campylobacter</i></b>	<b><i>Vibrio para.</i></b>	Cronobacter
Average read quality Q score for R1 and R2	>=30	>=30	>=30	>=30	>=30	>=30	>=30
Average coverage	>=30X	>=20X	>=40X	>=40X	>=20X	>=40X	>=20X
<i>De novo</i> assembly: Seq. length (Mbp)	~4.3-5.2	~2.7-3.2	~4.5-5.9	~4.0-5.0	~1.5-1.9	~4.8-5.5	~4-5
<i>De novo</i> assembly: no. contigs	<=300	<=300	<=500	<=650	<=300	<=300	<=500