

A



Nov 02, 2021

## © GCaMP6F Fluorescence in Ex Vivo Intestinal Preparation

Bryan Yoo<sup>1</sup>, Jessica Griffiths<sup>1</sup>, Sarkis Mazmanian<sup>1</sup>

<sup>1</sup>California Institute of Technology



dx.doi.org/10.17504/protocols.io.bzqap5se

## Mazmanian Lab

Jessica Griffiths

Protocol for GCaMP6F imaging in ex vivo intestinal tissue used in Yoo et al 2021

DOI

dx.doi.org/10.17504/protocols.io.bzqap5se

Bryan Yoo, Jessica Griffiths, Sarkis Mazmanian 2021. GCaMP6F Fluorescence in Ex Vivo Intestinal Preparation. **protocols.io** 

https://dx.doi.org/10.17504/protocols.io.bzqap5se

Center for Environmental and Microbial Interactions (CEMI) Grant ID: n/a	
Emerald Foundation Grant ID: n/a	
Heritage Medical Research Institute Grant ID: n/a	
Department of Defense Grant ID: PD160030	
National Institutes of Health Grant ID: GM007616 and DK078938	
Aligning Science Across Parkinson's Grant ID: ASAP-000375	
protocol ,	

Nov 02, 2021



1

## 54754

- 1 Small intestinal tissue was quickly harvested from ChAT-Cre mice, flushed and placed in oxygenated (95% O2, 5% CO2), ice cold Krebs-Henseleit solution for 1 hour followed by 15 min at room temperature.
- 2 A segment was cut along the mesenteric attachment, and pinned flat (mucosa facing down) on a Sylgard-lined recording chamber (Warner Instruments, PH1) in oxygenated Krebs-Henseleit solution.
- 3 C21 was added at 10nM and GCaMP6F fluorescence was detected on an upright microscope (Zeiss, Oberkochen, Germany-Examiner D1).