




AUG 02, 2023

Sending samples to the NCI Frederick EM Core

Bryce
Killingsworth¹

¹Translational Nanobiology Section, Laboratory of Pathology, Center for Cancer Research, National Cancer Institute, National Institutes of Health



Jennifer Jones

ABSTRACT

Protocol on how to ship EV sample to the NCI Frederick core.

OPEN  ACCESS



Protocol Citation: Bryce Killingsworth 2023. Sending samples to the NCI Frederick EM Core. **protocols.io** <https://protocols.io/view/sending-samples-to-the-nci-frederick-em-core-bjfdkji6>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working
We use this protocol and it's working

Created: Aug 06, 2020

Last Modified: Aug 02, 2023

PROTOCOL integer ID:
40133

Sending samples

- 1 If you are beginning a new project, you must put in a service request, select a CAN to draw funds from, and get approval from an AO, using the "New Request" button on the NCI Frederick EM Core's website:

<https://ncifrederick.cancer.gov/services/accessioning/services/labservices/area/5>

Pg 1: information about the request: fill out scope of work and number of samples

Pg2: select CCR investigator and tag emails of anyone you want to add as followers

Pg3: Daily Cruz-Abreu is the Fiscal Authority. Funding source is the current fiscal year Project ID number located on the "General Info and Links" tab of the lab OneNote page.

- 2 Ensure that Ziqui Wang or Kunio Nagashima at NCI Frederick is expecting your package delivery, and you have had communication about what type of images you expect.

- 3 The Frederick Courier no longer accepts requests via email or phone. You must use the Shipping Wizard to arrange pickup.

<https://ncifrederick.cancer.gov/cad/shippingwizard>

NCI Frederick Courier: 301-846-1175 (contact number only, do not call for pickup)

Core Address:

TO: Ziqui Wang

Electron Microscope Facility NCI-Frederick

EML at the ATRF

8560 Progress Drive


Wing D, Room D-2025

301-846-6064

- 4 Thaw the EV samples you intend to have imaged.
- 5 Transfer at least 50 μ L of each of your sample to 0.5 mL tubes for shipment. It is possible as little as 35 μ L could be imaged but any significant volume differences need to be discussed with the NCI core staff. Previously we have sent 100 μ L per sample.
- 6 Using a tube rack, other plastic supports, and tape, position the tubes in a plastic bag so that they are secure.

7 Place the tubes in a sytrofoam container with ice from the ice machine, not dry ice. Use a paper towel as a barrier between the ice and the samples so they don't partially freeze.

8 Create an item list using this template and print a copy to be shipped inside the box.



9 Tape up the box.

10 Affix a To and From label on the top of the box:

TO: Ziqiu Wang
Electron Microscope Facility NCI-Frederick
EML at the ATRF
8560 Progress Drive
Wing D, Room D-2025

FROM: Your Name
Jones Lab
Laboratory of Pathology
10 Center Dr
Bld 10, Rm B1B53
Bethesda MD,20892

11 You can wait for the courier to knock on the lab door, or place the package with a note for the courier outside the lab in the hallway.

12 Email Ziqui and Kunio to let them know the shipment is on the way.

nagashimak@mail.nih.gov

ziqui.wang@nih.gov