

NOV 12, 2023



## ( In vitro phosphorylation

Leonardo A Parra-Rivas<sup>1</sup>

<sup>1</sup>University of California, San Diego

Subhojit Roy Lab



Leonardo A Parra-Rivas University of California, San Diego

**ABSTRACT** 

In vitro phosphorylation





DOI:

dx.doi.org/10.17504/protocol s.io.ewov1q71ogr2/v1

**Protocol Citation:** Leonardo A Parra-Rivas 2023. In vitro phosphorylation. protocols.io https://dx.doi.org/10.17504/p rotocols.io.ewov1q71ogr2/v1

License: This is an open access protocol distributed under the terms of the Creative Commons Attribution License. which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status: Working** We use this protocol and it's working

Created: Nov 12, 2023

Last Modified: Nov 12,

2023

Oct 12 2023

## **PROTOCOL** integer ID:

90821

- 1 5  $\mu g$  of  $\alpha$ -syn GST fusion protein on beads were resuspended in a buffer containing 20 mM HEPES, 1 mM ATP, 10 mM MgCl<sub>2</sub>, 2 mM DTT, pH 7.4.
- The reaction was initiated by adding 1 ug of human recombinant PLK2 (Abcam, Cat#ab102108) and incubating at 30 °C with orbital shaking overnight (Eppendorf ThermoMixer F1.5).
- 3 Samples were then analyzed by NuPAGE and immunoblotted or alternatively analyzed by LC-MS/MS