



2 ▼

Apr 11, 2022

Suggested protocol for loading a DNA Ladder/marker V.2

[New England Biolabs¹](#)¹New England Biolabs

1

dx.doi.org/10.17504/protocols.io.zkxygx6zg8j2/v2**New England Biolabs (NEB)**Tech. support phone: **+1(800)632-7799** email: **info@neb.com****New England Biolabs**
New England Biolabs

This is the suggested protocol for use with:

Quick-Load[®] Purple 1 kb DNA Ladder ([N0552](#))
Quick-Load[®] Purple 100 bp DNA Ladder ([N0551](#))
Quick-Load[®] Purple 50 bp DNA Ladder ([N0556](#))
Quick-Load[®] 1 kb Extend DNA Ladder ([N3239](#))
Quick-Load[®] 1 kb DNA Ladder ([N0468](#))
Supercoiled DNA Ladder ([N0472S](#))
λ DNA-Mono Cut Mix ([N3019](#))
φX174 DNA-HaeIII Digest ([N3026](#))
pBR322 DNA-BstNI Digest ([N3031](#))
pBR322 DNA-MspI Digest ([N3032](#))
2-Log DNA Ladder (0.1-10.0 kb) ([N3200](#))
100 bp DNA Ladder ([N3231](#))
1 kb DNA Ladder ([N3232](#))
Low Molecular Weight DNA Ladder ([N3233](#))
50 bp DNA Ladder ([N3236](#))

DOI

dx.doi.org/10.17504/protocols.io.zkxygx6zg8j2/v2<https://www.neb.com/protocols/2018/08/03/suggested-loading-protocol-for-dna-ladders-and-markers>

New England Biolabs 2022. Suggested protocol for loading a DNA Ladder/marker. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.zkxygx6zg8j2/v2>
New England Biolabs



agarose gel, how to load a gel, loading a ladder, loading a marker


_____ protocol ,

May 10, 2020

Apr 11, 2022

36824

MATERIALS

 [Gel Loading Dye, Purple \(6X\), no SDS - 4.0 ml](#) **New England Biolabs Catalog #B7025S**

Please refer to the Safety Data Sheets (SDS) for health and environmental hazards.

This protocol is recommended for a 5 mm wide gel lane. The components of the mixture should be scaled up or down, depending on the width of the lane.

1



Prepare loading mixture (**6 µl** total volume):

Dilute only 1 µl of DNA Ladder at a time.

A	B
Distilled water (dH ₂ O)* or TE Buffer	4 µl
Gel Loading Dye, Purple (6X), no SDS	1 µl
DNA Ladder/Marker	1 µl
<i>Total Volume</i>	6 µl

*For multiple loads, dilution, and storage, use TE or other buffer of minimal ionic strength instead of water. DNA may denature if diluted and stored in dH₂O.



Mix gently by pipetting.



Load onto the agarose gel.