



2 ▼

Feb 20, 2022

Removal of Single-Stranded Extension using Mung Bean Nuclease (M0250) V.2

[New England Biolabs¹](#)¹New England Biolabs

1

dx.doi.org/10.17504/protocols.io.bb2riqd6**New England Biolabs (NEB)**Tech. support phone: **+1(800)632-7799** email: **info@neb.com****New England Biolabs**
New England Biolabs

This is the protocol for the removal of single-stranded extension using Mung Bean Nuclease (M0250).

DOI

dx.doi.org/10.17504/protocols.io.bb2riqd6<https://www.neb.com/protocols/0001/01/01/removal-of-single-stranded-extension>

New England Biolabs 2022. Removal of Single-Stranded Extension using Mung Bean Nuclease (M0250). **protocols.io**

<https://dx.doi.org/10.17504/protocols.io.bb2riqd6>[New England Biolabs](#)

single extension, removing single stranded extensions, mung bean nuclease



_____ protocol ,

Jan 31, 2020


Feb 20, 2022

32561

MATERIALS

[Mung Bean Nuclease - 1,500 units](#) **New England****Biolabs Catalog #M0250S**

Please see SDS (Safety Data Sheet) for hazards and safety warnings.

1 Suspend DNA ( **0.1 µg/µl**) in 1X *Mung Bean Nuclease Buffer* or 1X *NEBuffer 1.1, 2.1* or *CutSmart®*.

2 

Add 1.0 unit of *Mung Bean Nuclease* per µg DNA.

3 

Incubate at  **30 °C** for  **00:30:00** .

4 Inactivate the enzyme by phenol/chloroform extraction or by addition of SDS to 0.01%.

5 Recover the DNA by ethanol precipitation.