



VERSION 2
NOV 13, 2023

OPEN ACCESS



DOI:
dx.doi.org/10.17504/protocols.io.36wgq36b5lk5/v2

Protocol Citation: Willian Barela Costa 2023. Sea Water - Blue Treasure.
protocols.io
<https://dx.doi.org/10.17504/protocols.io.36wgq36b5lk5/v2>
Version created by Willian Barela Costa

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working
We use this protocol and it's working

🌐 Sea Water - Blue Treasure V.2

Willian Barela
Costa^{1,2}

¹FabLab Casa Firjan; ²@casafirjan



Willian Barela Costa

ABSTRACT

Pt 🇧🇷

Produção de água do mar para cultivo de microalgas utilizando a marca Blue Treasure.

25 - 35 Salinidade

En

Sea water production for microalgae cultivation using the Blue Treasure supplier.

25 - 35 Salinity

[M] 27 mg/mL

MATERIALS

Pt 🇧🇷

Sal do produto Blue Treasure
Água destilada

En



Salt from Blue Treasure product
Distilled water

Created: Nov 13, 2023



Last Modified: Nov 13, 2023

PROTOCOL integer ID:
90859

Metodologia - 1 litro

- 1 [M] 0.027 g/mL ->  Sample 27 g Sal do mar para  1 L água destilada
- 2 Agite até diluir
- 3 Certifique-se de que a salinidade esteja dentro da faixa desejada usando o refratômetro de salinidade. Neste protocolo, o foco é a salinidade 27.

Salting process - 1 liter

- 4 [M] 0.027 g/mL ->  Sample 27 g Sea Sal to  1 L Water distilled
- 5 Shake until diluted
- 6 Make sure the salinity is within the desired range using the salinity refractometer. In this protocol, the focus is 27 salinity.