

NOV 02, 2023

# FTIR

## Alexander Röntgen<sup>1</sup>

<sup>1</sup>University of Cambridge



Alexander Röntgen

### **ABSTRACT**

Protein FTIR protocol.

### **MATERIALS**

Vertex 70 FTIR spectrometer with a  $\it DiamondATR$  unit and a deuterated lanthanum  $\it \alpha$ -alanine-doped triglycine sulphate detector (Bruker)





#### DOI:

dx.doi.org/10.17504/protocol s.io.dm6gp3xqpvzp/v1

**Protocol Citation:** Alexande r Röntgen 2023. FTIR. **protocols.io** https://dx.doi.org/10.17504/protocols.io.dm6gp3xqpvzp/v1

License: This is an open access protocol distributed under the terms of the Creative Commons
Attribution License, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working We use this protocol and it's working

Created: Oct 31, 2023

Last Modified: Nov 02,

2023

# **Background scan**

- 1 Spot 5 μL of buffer on prism
  - Acquire background scan

## Sample scan

- 2 Spot 5 μL of protein solution on prism
  - Acquire sample scan