

Apr 18, 2024 Version 2

Glucose Tolerance Test V.2

DOI

dx.doi.org/10.17504/protocols.io.261geddodv47/v2

Sabina Marciano¹, Roberta Marongiu¹

¹Weill Cornell Medicine

ASAP Collaborative Rese...



Eileen

Eileen Ruth Torres

Weill Cornell Medicine

OPEN  ACCESS



DOI: **dx.doi.org/10.17504/protocols.io.261geddodv47/v2**

Protocol Citation: Sabina Marciano, Roberta Marongiu 2024. Glucose Tolerance Test. **protocols.io**

<https://dx.doi.org/10.17504/protocols.io.261geddodv47/v2>Version created by **[Eileen Ruth Torres](#)**

License: This is an open access protocol distributed under the terms of the **[Creative Commons Attribution License](#)**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: April 18, 2024

Last Modified: April 18, 2024

Protocol Integer ID: 98408

Keywords: ASAPCRN

Funders Acknowledgement:

Aligning Science Across

Parkinson's


Grant ID: 020608



Abstract

The glucose tolerance test is performed to determine how quickly the glucose is cleared from the blood. It is used to test for **diabetes** or **insulin resistance**.



- 1 Single cage the mice for 1 week in advance
- 2 Weight the mice
- 3 Fast mice for  06:00:00 6h
- 4 After fasting, measure glucose with a glucose meter
- 5 Perform an intraperitoneal injection of 2g / kg body weight of glucose (20% D-glucose stock solution dissolving 2g of glucose in 10ml saline and give 10ul per gram of body weight).
- 6 Collect blood samples for glucose measurement at 15, 30, 45, 60, 90 and 120 min post glucose injection