

Jun 23, 2024

Synapse Quantification with Puncta Analyzer

 In 1 collection

DOI

dx.doi.org/10.17504/protocols.io.8epv5rnddg1b/v1

Justin T Savage¹

¹Duke University

ASAP Collaborative Rese...

Eroglu_Lab



Justin T Savage

Duke University

OPEN  ACCESS



DOI: **dx.doi.org/10.17504/protocols.io.8epv5rnddg1b/v1**

Protocol Citation: Justin T Savage 2024. Synapse Quantification with Puncta Analyzer. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.8epv5rnddg1b/v1>

License: This is an open access protocol distributed under the terms of the **[Creative Commons Attribution License](#)**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

Created: June 23, 2024

Last Modified: June 23, 2024

Protocol Integer ID: 102277

Keywords: ASAPCRN

Abstract

Protocol for installing Puncta Analyzer and using it to analyze synapse images.



Install Puncta Analyzer

- 1 Go to <https://github.com/toddstavish/puncta-analyzer> and follow the following instructions for installing Puncta Analyzer
- 2 Install JDK
- 3 Install maven
- 4 Fork the puncta analyzer repo
- 5 Execute 'mvn compile' in project directory (where the pom.xml is located)
- 6 Create a 'PunctaAnalyzer' directory in the your ImageJ plugins directory
- 7 Copy all of the class files (bytecode) in target/classes to plugins/PunctaAnalyzer

Running Puncta Analyzer

- 8 Z-project images if necessary
- 9 Open an image in FIJI and use Ctrl+A to select all or select a specific region of interest
- 10 Run Puncta Analyzer using Analyze>Puncta Analyzer
- 11 Select Red channel and Green channel boxes to analyze those channels. Also select background subtraction for the red and green channels. Click set results file and save results and select a folder to save to.



- 12 Click OK with the default 50 pixel rolling ball radius for subtract background
- 13 Threshold the red channel using the slider
- 14 Click OK with the default 50 pixel rolling ball radius for subtract background
- 15 Threshold the red channel using the slider
- 16 Enter any desired restrictions on the analyze particles function such as a minimum pixel size.
- 17 Click yes to save result to file
- 18 Close all windows and proceed to the next image