



APR 01, 2024

🌐 CODA: shorthand for calling functions | HuBMAP | JHU-TMC

Kyu Sang Han¹, Pei-Hsun Wu¹, Joel Sunshine², Ashley Kiemen²,
Sashank Reddy², Denis Wirtz^{1,2}

¹Johns Hopkins University; ²Johns Hopkins Medicine

Human BioMolecular Atlas Program (HuBMAP) Method Development Community

TMC - Johns Hopkins University

OPEN  ACCESS



Kyu Sang Han
Johns Hopkins University



DOI:

dx.doi.org/10.17504/protocols.io.5qpvoknkb14o/v1

Protocol Citation: Kyu Sang Han, Pei-Hsun Wu, Joel Sunshine, Ashley Kiemen, Sashank Reddy, Denis Wirtz 2024. CODA: shorthand for calling functions | HuBMAP | JHU-TMC. [protocols.io](https://dx.doi.org/10.17504/protocols.io.5qpvoknkb14o/v1) <https://dx.doi.org/10.17504/protocols.io.5qpvoknkb14o/v1>

MANUSCRIPT CITATION:

Kiemen, A.L., Braxton, A.M., Grahn, M.P. *et al.* CODA: quantitative 3D reconstruction of large tissues at cellular resolution. *Nat Methods* **19**, 1490–1499 (2022). <https://doi.org/10.1038/s41592-022-01650-9>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working
We use this protocol and it's working

Created: Mar 28, 2024

Last Modified: Apr 01, 2024

PROTOCOL integer ID: 97516

Keywords: CODA, machinelearning, deeplearning

Funders Acknowledgement:

Institute of Arthritis and Musculoskeletal and Skin Diseases
Grant ID: U54AR081774
National Cancer Institute
Grant ID: U54CA143868

ABSTRACT

To downsample ndpi or svf images to 10x, 5x, and 1x tifs, use this function:

create_downsampled_tif_images

or try Openslide in python

To calculate registration on the low resolution (1x) images

1. calculate the tissue area and background pixels using this function:

calculate_tissue_ws

2. calculate the registration transforms:

calculate_image_registration

To build a 3D tissue volume using semantic segmentation:

1. generate manual annotations in Aperio imagescope

2. apply the deep learning function to train a model and segment the high resolution (5x or 10x) images:

train_image_segmentation

To apply the registration to segmented images:

apply_image_registration

To build a 3D tissue matrix from registered, classified images:

build_tissue_volume

To build a 3D cell volume containing nuclear coordinates:

1. Build a mosaic image containing regions of many whole-slide images for cell detection optimization:

make_cell_detection_mosaic

2. Manually annotate the mosaic image to get the 'ground-truth' number of cell nuclei:

manual_cell_count

3. Determine cell detection parameters using the manual annotations on the mosaic image:

get_nuclear_detection_parameters

4. Deconvolve the high-resolution (5x or 10x) H&E images before applying the cell detection algorithm:

deconvolve_histological_images

5. Detect cells on the hematoxylin channel of the high-resolution images:

cell_detection

6. Apply the registration to the cell coordinates:

register_cell_coordinates

7. Build a 3D cell coordinate matrix corresponding to the 3D tissue matrix:

build_cell_volume