

DEC 01, 2023

OPEN BACCESS



DOI:

dx.doi.org/10.17504/protocol s.io.bp2l6xybzlqe/v1

Protocol Citation: Nandie Elhadidy 2023. Nephele DADA2 Workflow. protocols.io https://dx.doi.org/10.17504/p rotocols.io.bp2l6xybzlqe/v1

License: This is an open access protocol distributed under the terms of the Creative Commons Attribution License, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working We use this protocol and it's working

Created: Dec 01, 2023

Last Modified: Dec 01,

2023

Nephele DADA2 Workflow

Nandie Elhadidy¹

¹North Carolina State University



Nandie Elhadidy

ABSTRACT

A protocol outlining the workflow of utilizing DADA2 on Nephele with a provided practice dataset. DADA2 is a pipeline that provides species assignments using Amplicon Sequence Variants or ASVs which compare variants to a predetermined sequence through assumption that the most amplified sequence is the correct one.

PROTOCOL integer ID:

91672

Downloading the Dataset

- 1 To download the practice dataset provided by Nephele, visit: https://nephele.niaid.nih.gov/user_guide_data/
- 2 Download the first set of sequences listed as well as the mapping file.
- 2.1 The reference provided for this set of sequences is listed as: Experimental Microbial Dysbiosis Does Not Promote Disease Progression in SIV-Infected Macaques. NCBI BioProject: PRJNA417022
- 3 Once downloaded, unzip the downloaded file and ensure all .fastq.gz files are present.

Running DADA2

- 4 Click the **Home** tab.
- 5 Choose the **Analyze** option.
- 6 Under **Amplicon Metagenomics**, select "Start job" option under **DADA2**.

7 Select **Paired End FASTQ**, as this is a paired-end dataset. 8 Click "Upload from my computer" and add all .fastq.gz files from the downloaded dataset. Start upload 9 Select "Next" and upload the corresponding Mapping File. 10 After reaching the Pipeline Selection page, click **Select** and enter a description of the job. 11 Keep Filter and Trim, Merge Pairs, and Analysis settings defaulted. Do not check Ion Torrent Data. Click Validate and Submit *Note: This is for following the exact tutorial instructions. If you would like to make changes to these pre-sets, proceed with these changes. 12 Nephele will send an email upon job completion. **Analysis of the Output** 13 From the My Workspace tab, hover to the view results icon under Actions for the completed DADA2 job. 14 Under Reported Statistics, Reported Files, Reported Images, Reported Tables, and Reported Charts and explore outputs.

14.1 Reported files of note include: Rooted tree, Rarefaction curve, Alpha diversity, PCoA binomial, BCoA Bray-Curtis, and ASV heatmap.

14.2 Reported images of note include: Quality profiles of forward and reverse reads

14.3 Reported tables of note include: Read counts by sample and Top species table

14.4 Reported charts of note include: Reads in vs. reads out and Track reads