



Version 2 ▾

Jun 16, 2021

# Human Kidney, Urinary Tract, and Lung Cell Type Mapping Pipeline for the Human Biomolecular Atlas Program (HuBMAP) V.2

Blue Lake<sup>1</sup>, Elizabeth Duong<sup>1</sup>, Dinh H H. Diep<sup>1</sup>, Xin Sun<sup>1</sup>, Peter K<sup>2</sup>, Jim Hagood<sup>3</sup>, Gloria S Pryhuber<sup>4</sup>, Sanjay Jain<sup>5</sup>, Kun Zhang [University Of California]<sup>1</sup>

<sup>1</sup>University of California, San Diego; <sup>2</sup>Harvard Medical School; <sup>3</sup>University of North Carolina at Chapel Hill;

<sup>4</sup>University of Rochester Medical Center; <sup>5</sup>Washington University, Saint Louis

1 Works for me



Share

[dx.doi.org/10.17504/protocols.io.bvvin64e](https://dx.doi.org/10.17504/protocols.io.bvvin64e)

Human BioMolecular Atlas Program (HuBMAP) Method Development Community  
Tech. support email: [Jeff.spraggins@vanderbilt.edu](mailto:Jeff.spraggins@vanderbilt.edu)

Blue Lake  
UCSD

## ABSTRACT

This document provides an overview of the protocols used by the Human Kidney, Urinary Tract, and Lung Mapping Center (KULMAP) for generation of spatially resolved multomic cell type organ maps within HuBMAP. This involves sequencing of the transcriptomes and epigenomes of dissociated single cells in a massively parallel manner, which will then inform on a highly multiplexed RNA *in situ* imaging method (DART-FISH) for spatial mapping of hundreds of molecular targets in the tissue sections, at a subcellular resolution. These spatial molecular maps will serve as scaffolds for computational registration of cell types and the associated transcriptome/chromatin maps to the tissue space.

## DOI

[dx.doi.org/10.17504/protocols.io.bvvin64e](https://dx.doi.org/10.17504/protocols.io.bvvin64e)

## PROTOCOL CITATION

Blue Lake, Elizabeth Duong, Dinh H H. Diep, Xin Sun, Peter K, Jim Hagood, Gloria S Pryhuber, Sanjay Jain, Kun Zhang [University Of California 2021. Human Kidney, Urinary Tract, and Lung Cell Type Mapping Pipeline for the Human Biomolecular Atlas Program (HuBMAP). **protocols.io**  
<https://dx.doi.org/10.17504/protocols.io.bvvin64e>  
Version created by Blue Lake

## LICENSE

———— This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

## CREATED

Jun 16, 2021

## LAST MODIFIED

Jun 16, 2021

## PROTOCOL INTEGER ID

50826

- 1 Collection and processing of human tissues.

Bladder, Ureter, Kidney: [dx.doi.org/10.17504/protocols.io.568g9hw](https://doi.org/10.17504/protocols.io.568g9hw)  
Lung airways and parenchyma: [dx.doi.org/10.17504/protocols.io.bjuxknxn](https://doi.org/10.17504/protocols.io.bjuxknxn) and  
[dx.doi.org/10.17504/protocols.io.bjtnknme](https://doi.org/10.17504/protocols.io.bjtnknme)

- 2 Isolation of single nuclei for omic assays.  
Bladder, Ureter, Kidney: [dx.doi.org/10.17504/protocols.io.ufketkw](https://doi.org/10.17504/protocols.io.ufketkw)  
Lung airways and parenchyma: [dx.doi.org/10.17504/protocols.io.bh26j8he](https://doi.org/10.17504/protocols.io.bh26j8he)
- 3 Running of multimodal single cell/nucleus assays.  
10X snRNA-seq: [dx.doi.org/10.17504/protocols.io.86khzcw](https://doi.org/10.17504/protocols.io.86khzcw)  
10X snATAC-seq: [dx.doi.org/10.17504/protocols.io.bvssn6ee](https://doi.org/10.17504/protocols.io.bvssn6ee)  
SNARE-Seq2: [dx.doi.org/10.17504/protocols.io.be5gig3w](https://doi.org/10.17504/protocols.io.be5gig3w)
- 4 Spatial in situ mapping of multiomic profiles.  
DART-FISH: Protocols coming soon.