



Jul 31, 2020

NEBNext End Repair Module E6050

Isabel Gautreau¹¹New England Biolabs
1 Works for me dx.doi.org/10.17504/protocols.io.k64czgw

New England Biolabs (NEB)

Tech. support phone: +1(800)632-7799 email: info@neb.com
 Isabel Gautreau
New England Biolabs

EXTERNAL LINK

<https://www.neb.com/-/media/catalog/datacards-or-manuals/manuale6050.pdf>

DOI

dx.doi.org/10.17504/protocols.io.k64czgw

PROTOCOL CITATION

Isabel Gautreau 2020. NEBNext End Repair Module E6050. **protocols.io**
dx.doi.org/10.17504/protocols.io.k64czgw



EXTERNAL LINK

<https://www.neb.com/-/media/catalog/datacards-or-manuals/manuale6050.pdf>

LICENSE

This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

Dec 07, 2017

LAST MODIFIED

Jul 31, 2020

PROTOCOL INTEGER ID

9148

STEPS MATERIALS

NAME	CATALOG #	VENDOR
NEBNext End Repair Enzyme Mix	E6051	New England Biolabs
NEBNext End Repair Reaction Buffer:	E6052	New England Biolabs

Mix the following components in a sterile microfuge tube:

- Fragmented DNA: variable
NEBNext End Repair Reaction Buffer (10X): 10 µl
NEBNext End Repair Enzyme Mix: 5 µl
Sterile H₂O for a final volume of 100 µl: variable

Total volume: 100 µl



NEBNext End Repair Enzyme Mix

by New England Biolabs

Catalog #: E6051



NEBNext End Repair Reaction Buffer:

by New England Biolabs

Catalog #: E6052

- 2 Incubate in a thermal cycler for 30 minutes at 20°C.
- 3 Purify DNA sample on one column.