




AUG 14, 2023

Protocol details for autopsy specimens (WUSTL)

Reyka Jayasinghe¹

¹Washington University in St. Louis



Reyka Jayasinghe

ABSTRACT

protocol details for autopsy specimens

OPEN  ACCESS



DOI:

dx.doi.org/10.17504/protocols.io.e6nvwdd92lmk/v1

Protocol Citation: Reyka Jayasinghe 2023. Protocol details for autopsy specimens (WUSTL). **protocols.io** <https://dx.doi.org/10.17504/protocols.io.e6nvwdd92lmk/v1>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working
We use this protocol and it's working

Created: Aug 14, 2023

Last Modified: Aug 14, 2023

Sample Collection Strategy

- 1 The tissue is transported on ice to the lab for processing and storage.

Larger samples are sectioned into small segments and placed into 2 mL cryovial before dropping into liquid nitrogen.

Vials are immediately placed in -80° freezer