



Dec 11, 2020

m9_media

Elizabeth Fozo¹¹In-house protocol**1***Works for me*

This protocol is published without a DOI.

**Eadewunm**

ABSTRACT

M9 media protocols

THIS PROTOCOL ACCOMPANIES THE FOLLOWING PUBLICATION

[Cold Spring Harbor Protocols](#)

PROTOCOL CITATION

Elizabeth Fozo 2020. m9_media. **protocols.io**
<https://protocols.io/view/m9-media-bqmtmu6n>

MANUSCRIPT CITATION please remember to cite the following publication along with this protocol

[Cold Spring Harbor Protocols](#)

KEYWORDS

Cold Spring Harbor Protocols, M9 media

LICENSE

————— This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

Dec 11, 2020

LAST MODIFIED

Dec 11, 2020

PROTOCOL INTEGER ID

45459

GUIDELINES

M9 media

Cold Spring Harbor Protocols

DISCLAIMER:

DISCLAIMER: THIS WORK IS IN PROGRESS. IT IS FOR INFORMATIONAL PURPOSES ONLY; USE AT YOUR OWN RISK The protocol content here is for informational purposes only and does not constitute legal, medical, clinical, or safety advice, or otherwise; content added to [protocols.io](#) is not peer-reviewed and may not have undergone a formal approval of any kind. Information presented in this protocol should not substitute for independent professional judgment, advice, diagnosis, or treatment. Any action you take or refrain from taking using or relying upon the information presented here is strictly at your own risk. You agree that neither the Company nor any of the authors, contributors, administrators, or anyone else associated with [protocols.io](#), can be held responsible for your use of the information contained in or linked to this protocol or any of our Sites/Apps and Services.

ABSTRACT

M9 media protocols

M9 minimal media: low osmolarity media for E. coli, resulting in slower growth rate

- 1 Mix the ingredients (forms white precipitate but is dissolved after steering), filter sterilize, and store at room temperature.

Composition	Required	Stock	Vol(ml)/L	Vol(ml)/500mL
M9 Salts	1X	5X	200	100
Magnesium sulfate	2mM	1M	2	1
Calcium chloride	0.1mM	1M	0.1	0.05
Vitamin B1	1ug/mL	10mg.mL	0.1	0.05
Carbon source			as required	as required
Distilled water			up to 1000mL	up to 500mL

Cold Spring Harbor Protocols:

- 2 http://cshprotocols.cshlp.org/content/2006/1/pdb.rec8146.full?text_only=true