




APR 15, 2024

In vivo Ephys: Spike Sorting Protocol

 In 2 collections

Sasha Burwell¹

¹Duke University

ASAP Collaborative Research Network



Sasha Burwell
Duke University

OPEN  ACCESS



ABSTRACT

This protocol details the spike sorting procedures used on the collected electrode recording data.

DOI:

dx.doi.org/10.17504/protocols.io.8epv5r1b4g1b/v1

Protocol Citation: Sasha Burwell 2024. In vivo Ephys: Spike Sorting Protocol. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.8epv5r1b4g1b/v1>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

Created: Apr 02, 2024

Last Modified: Apr 15, 2024

Spike Sorting Protocol

1 For more information on setting up and using Spyking Circus, see: <https://spyking-circus.readthedocs.io/en/latest/code/index.html> .

1.1 Make a copy of the continuous.dat data file from your recording into a new “sorted_data” folder.

Note

Important note: Make sure to place a **copy** the continuous.dat file in this folder, and not the original, as the file will change as a result of spike sorting and you want to preserve the original raw data. This means you will also need to make a new copy the raw continuous.dat file if you need to redo the spike sorting analysis for any reason.

2 Add a .params file with the same name as the .dat file to the folder (i.e., continuous.params). The continuous.params used for these experiments:



```
[data]
file_format      = raw_binary
data_offset      = 0
mapping          = path\open_ephys_16prb.prb      # Mapping of the electrode

#see prb file below
suffix           =                # Suffix to add to generated files
data_dtype       = int16          # Type of the data
dtype_offset     = auto           # Padding for data (if auto: int16 is 0)
skip_artefact    = False         # Skip part of signals with large fluctuations
sampling_rate    = 30000         # Sampling rate of the data [Hz]
stationary       = True          # Should be False for long recordings
radius           = auto          # Radius [in um] (if auto, read from the prb
file)
alignment        = True          # Realign the waveforms by oversampling
global_tmp       = True
multi-files      = False         # If several files mydata_0,1,...,n.dat should be
processed together

[detection]
spike_thresh     = 4
N_t              = 2 or 3        # Width of the templates [in ms]. We used both 2
and 3 in different sortings to extract all the templates included.
peaks            = positive      # can be positive, negative, or both

[filtering]
cut_off          = 250           # Cut off freq for the butterworth filter [Hz]
filter           = True         # If True, then a low-pass filtering is performed

[whitening]
chunk_size       = 60           # Size of the data chunks [in s]
safety_time      = 1           # Temporal zone around which templates are
isolated [in ms]
temporal         = False        # Perform temporal whitening
spatial          = True         # Perform spatial whitening
max_elts         = 10000        # Max number of events per electrode
nb_elts          = 0.8          # Fraction of max_elts that should be obtained
per electrode [0-1]
output_dim       = 5            # Can be in percent of variance explain, or num
of dimensions for PCA on waveforms

[clustering]
extraction       = median-row    # Can be either median-row (default), median-
pca, mean-pca, mean-row, or quadratic
```

```

safety_space    = True          # If True, we exclude spikes in the vicinity of a
selected spikes
safety_time     = auto          # Temporal zone around which templates are
isolated [in ms]
max_elts        = 15000         # Max number of events per electrode (should be
compatible with nb_elts)
nb_elts         = 0.8           # Fraction of max_elts that should be obtained
per electrode [0-1]
nclus_min       = 0.0075        # Min number of elements in a cluster (given in
percentage)
max_clusters    = 20            # Maximal number of clusters for every electrodes
nb_repeats      = 3             # Number of passes used for the clustering
smart_search    = 0             # Parameter for the smart search [0-1]. The
higher, the more strict
sim_same_elec   = 0.75          # Distance within clusters under which they are
re-merged
cc_merge        = 1             # If CC between two templates is higher, they are
merged
noise_thr       = 0.8           # Minimal amplitudes are such than
amp*min(templates) < noise_thr*threshold
make_plots      = png           # Generate sanity plots of the clustering
remove_mixture  = True          # At the end of the clustering, we remove
mixtures
of templates

[fitting]
chunk           = 1             # Size of chunks used during fitting [in second]
gpu_only        = False         # Use GPU for computation of b's AND fitting
amp_limits      = (0.3, 5)      # Amplitudes for the templates during spike
detection
amp_auto        = True          # True if amplitudes are adjusted automatically
for every templates
refractory      = 0             # Refractory period, in ms [0 is None]
max_chunk       = inf           # Fit only up to max_chunk

[merging]
cc_overlap      = 0.5           # Only templates with CC higher than cc_overlap
may be merged
cc_bin          = 2             # Bin size for computing CC [in ms]

[extracting]
safety_time     = 3             # Temporal zone around which spikes are isolated
[in ms]
max_elts        = 1000          # Max number of collected events per templates
output_dim      = 5             # Percentage of variance explained while

```

```
performing PCA
cc_merge      = 1          # If CC between two templates is higher, they are
merged
noise_thr     = 0.8        # Minimal amplitudes are such than
amp*min(templates) < noise_thr*threshold

[noedits]
filter_done   = False      # Automatically edited - make sure this is False
before running the sorter
```

3

The .params file needs to have correct path to a .prb file (the mapping of the electrode you are using). We recorded with bundles of 16 individual electrodes, and used the following .prb file:



```
total_nb_channels = 16
radius            = 80
channel_groups = {
  1: {
    'channels': list(range(16)),
    'graph' : [],
    'geometry': {
      0: (0, 300),
      1: (100, 300),
      2: (200, 300),
      3: (300, 300),
      4: (0, 200),
      5: (100, 200),
      6: (200, 200),
      7: (300, 200),
      8: (0, 100),
      9: (100, 100),
      10: (200, 100),
      11: (300, 100),
      12: (0, 0),
      13: (100, 0),
      14: (200, 0),
      15: (300, 0),
    }
  }
}
```

- 4 Launch the Spyking Circus GUI from the anaconda terminal.
- 5 Run the sorting in the Main Algorithm tab.
 - 5.1 Browse → select continuous.dat file.
 - 5.2 Increase CPUs as desired, then click “Run”.
- 6 When the next window pops up, click “Suggest pairs”, “Select pairs”, then “Merge and finalize”.
- 7 To view results: View Results tab → Browse → select the analyzed continuous.dat file.
 - 7.1 Options:
 1. Matlab GUI
 2. Extension = “merged”
 - 7.2 Click “Run”.

7.3 See <https://spying-circus.readthedocs.io/en/latest/GUI/sorting.html> for detailed documentation on what each panel of the GUI does.

8 Sorting:

8.1 Set bin size to 1, click “Big ISI”, and update the scale so you can see whole waveforms (zooming out usually works).

8.2 Go through all templates.

8.3 First pass: Kill a cell if it has above a 1% incidence of refractory period violations.

Note

RPV: % of spike pairs where your assumption about the refractory period are violated.

8.4 Second pass: merge templates that you think are the same cell, using the two template view.

- Criteria used: similar waveforms, cross correlation dropping to zero at zero, coefficient of similarity close to 1.
- Sometimes you have to merge more than two templates for one cell.
- Merge template with fewer spikes into the one with more spikes.

8.5 Third pass: kill units that don’t meet requirements (too sparse, not spiking consistently across the whole session, auto correlation doesn’t go to zero, RPV has gone above 1%).

- 9 Add suffix (ex, 'v1' for first attempt at sorting) and save.