


Jul 24, 2020

## SOLUTION- 01 - Sodium Chloride (NaCl) solution

Marco Cosentino<sup>1</sup>, Elisa Storelli<sup>1</sup>, Alessandra Luini<sup>1</sup>, Massimiliano LM Legnaro<sup>1</sup>, Emanuela Rasini<sup>1</sup>, Marco Ferrari<sup>1</sup>, Franca Marino<sup>1</sup>

<sup>1</sup>Center for Research in Medical Pharmacology, University of Insubria (Varese, Italy)

1 Works for me dx.doi.org/10.17504/protocols.io.biyvkfw6



DOI

dx.doi.org/10.17504/protocols.io.biyvkfw6

DOCUMENT CITATION

Marco Cosentino, Elisa Storelli, Alessandra Luini, Massimiliano LM Legnaro, Emanuela Rasini, Marco Ferrari, Franca Marino 2020. SOLUTION-01 - Sodium Chloride (NaCl) solution. **protocols.io** dx.doi.org/10.17504/protocols.io.biyvkfw6

LICENSE

This is an open access document distributed under the terms of the Creative Commons Attribution License, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

Jul 24, 2020

LAST MODIFIED

Jul 27, 2020

DOCUMENT INTEGER ID

39669

NaCl composition g/I (M):

NaCl 8.770 (0.15)

Add ultrapure water to 1 liter.

Storage: § 4 °C Fridge 1 - (Room TSO8).

NaCl code: S9625, Sigma

protocols.io
1
07/24/2020

Citation: Marco Cosentino, Elisa Storelli, Alessandra Luini, Massimiliano LM Legnaro, Emanuela Rasini, Marco Ferrari, Franca Marino (07/24/2020). SOLUTION- 01 - Sodium Chloride (NaCl) solution. <a href="https://dx.doi.org/10.17504/protocols.io.biyvkfw6">https://dx.doi.org/10.17504/protocols.io.biyvkfw6</a>