




Aug 03, 2022

Processing of Radial K-space DW-MRI Data

Miguelrj¹, Hee Kwon Song¹, Stephen Pickup¹, Mamtaaryagupta¹,
Rong Zhou¹

¹University of Pennsylvania, Department of Radiology

1 Works for me



Share

dx.doi.org/10.17504/protocols.io.j8nlkkwj1l5r/v1



Miguelrj

ABSTRACT

The protocol includes reconstruction of diffusion weighted images from radial k-space data and using a 3-parameter fit to derive metrics of ADC and kurtosis index.

DOI

dx.doi.org/10.17504/protocols.io.j8nlkkwj1l5r/v1

PROTOCOL CITATION

Miguelrj , Hee Kwon Song, Stephen Pickup, Mamtaaryagupta , Rong Zhou 2022.
Processing of Radial K-space DW-MRI Data. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.j8nlkkwj1l5r/v1>



FUNDERS ACKNOWLEDGEMENT

National Cancer Institute, NIH
Grant ID: U24CA231858

LICENSE

————— This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

Feb 23, 2022

LAST MODIFIED

Aug 03, 2022

1 Image reconstruction

Radially acquired diffusion-weighted images (DWIs) are reconstructed using the following steps:

1. Apply zero-order phase correction to each radial spoke using the average phase offset at the center slice of the lowest b-value image
2. Zerofill k-space by a factor of 2 to double field of view
3. Multiply signal of each point by its respective area on a Voronoi diagram of the points (including added zerofill points) in k-space
4. Re-grid each radially defined point to its nearest Cartesian coordinate using its Kaiser-Bessel index
5. Apply Fourier transform to now Cartesian-defined k-space

This process can be easily conducted using the **CIRP Python library** [link]. Once the library is installed, the DWI reconstruction class can be imported using

```
from CIRP.DWI_processing import DWI_reconstructor
```

To instantiate the class, the image dimensions, b-values, readout resolution, and number of views must be defined:

```
img_size = [16, 96, 96] # [slices, yres, xres]
bvalues = [10, 535, 1070, 1479, 2141] # b-values
xres_ro = 128 # readout points per view
n_views = 403

reconstructor = DWI_reconstructor(image_size=img_size,
bvalues=bvalues, xres_ro=xres_ro, views = nviews)
```