



OCT 12, 2023

OPEN ACCESS



DOI:

[dx.doi.org/10.17504/protocols.io.81wgbxxpnlpk/v1](https://dx.doi.org/10.17504/protocols.io.81wgbxxpnlpk/v1)

**Protocol Citation:** Sabina Marciano, Roberta Marongiu 2023. Cylinder Test.

**protocols.io**

<https://dx.doi.org/10.17504/protocols.io.81wgbxxpnlpk/v1>

**License:** This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working  
We use this protocol and it's working

**Created:** Aug 14, 2023

**Last Modified:** Oct 12, 2023


## Cylinder Test

Sabina

Marciano<sup>1</sup>, Roberta Marongiu<sup>1</sup>

<sup>1</sup>Weill Cornell Medicine

ASAP Collaborative Research Network



Dylan J Barbera

Weill Cornell Medicine


### ABSTRACT


The Cylinder test is used to evaluate locomotor asymmetry in rodent models

### MATERIALS

- Glass cylinder 20 cm in diameter 30 cm high
- Angled mirror
- 70% ethanol for cleaning between each session

**Keywords:** ASAPCRN

- 1 Place the mouse in the transparent cylinder for  00:05:00 5m
- 2 Place an angled mirror at the back of the cylinder to allow observation of behavior from every angle

- 3 Record mouse behavior for  00:05:00 5m

Main outcome measures: Left/Right contacts, Number of rears, L/R hind-paw contacts, Total time rearing, Time on the ground