



2 ▼

Oct 29, 2021

# A simple and economic protocol for efficient *in vitro* fertilization using cryopreserved mouse sperm V.2

PLOS One ✓

Magdalena Wigger<sup>1,2</sup>, Simon E Tröder<sup>1,2</sup>, Branko Zevnik<sup>1,2</sup><sup>1</sup>Cluster of Excellence Cellular Stress Responses in Aging-associated Diseases (CECAD), Faculty of Medicine and University Hospital Cologne, University of Cologne, Cologne, Germany;<sup>2</sup>in vivo Research Facility, Faculty of Medicine and University Hospital Cologne, University of Cologne, Cologne, Germany

1



dx.doi.org/10.17504/protocols.io.bx2spqee

PLOS ONE Lab Protocols

Tech. support email: [plosone@plos.org](mailto:plosone@plos.org)

Magdalena Wigger

The advent of genome editing tools like CRISPR/Cas has substantially increased the number of genetically engineered mouse models in recent years. In support of refinement and reduction, sperm cryopreservation is advantageous compared to embryo freezing for archiving and distribution of such mouse models. The *in vitro* fertilization using cryopreserved sperm from the most widely used C57BL/6 strain has become highly efficient in recent years due to several improvements of the procedure. However, purchase of the necessary media for routine application of the current protocol poses a constant burden on budgetary constraints. In-house media preparation, instead, is complex and requires quality control of each batch. Here, we describe a cost-effective and easily adaptable approach for *in vitro* fertilization using cryopreserved C57BL/6 sperm. This is mainly achieved by modification of an affordable commercial fertilization medium and a step-by-step description of all other necessary reagents. This protocol is compatible with frozen sperm from all major repositories and the IVF can easily be adapted to accommodate freshly harvested sperm.

Protocol scheme.png

DOI

[dx.doi.org/10.17504/protocols.io.bx2spqee](https://dx.doi.org/10.17504/protocols.io.bx2spqee)<https://journals.plos.org/plosone/article?id=10.1371/journal.pone.0259202>

Magdalena Wigger, Simon E Tröder, Branko Zevnik 2021. A simple and economic protocol for efficient *in vitro* fertilization using cryopreserved mouse sperm. **protocols.io**

<https://dx.doi.org/10.17504/protocols.io.bx2spqee>

Magdalena Wigger



protocol

Wigger M, Tröder SE, Zevnik B (2021) A simple and economic protocol for efficient *in vitro* fertilization using cryopreserved mouse sperm. PLOS ONE 16(10): e0259202.

<https://doi.org/10.1371/journal.pone.0259202>Sperm, Cryopreservation, Embryos, Genetics, Genome editing, *In vitro* fertilization, Mouse model

protocol ,

Sep 07, 2021





- Prepare a 100x  $\text{CaCl}_2 \times 2\text{H}_2\text{O}$  stock solution (310 mM) by dissolving 0.4558 g of  $\text{CaCl}_2 \times 2\text{H}_2\text{O}$  in 10 ml of embryo-grade water
- Filter the solution through a 0.22  $\mu\text{m}$  filter and store aliquots at  $-20^\circ\text{C}$  for a maximum of 6 months
- On the day of IVF thaw an aliquot of  $\text{CaCl}_2$  at room temperature
- Add 150  $\mu\text{l}$  of 100x  $\text{CaCl}_2$  to 15 ml of RVF medium and mix gently
- Place 1 ml of RVF medium supplemented with  $\text{CaCl}_2$  in a tube containing 30.7 mg of GSH and vortex
- Add 50  $\mu\text{l}$  of this solution to 5 ml of RVF medium supplemented with  $\text{CaCl}_2$ , mix gently and filter using 0.22  $\mu\text{m}$  syringe end filter in order to obtain mRVF medium

Note: HTF medium (e.g. Merck MR-070-D) may likely be used instead of RVF and supplemented as described above.

Use **M2** medium for handling of embryos outside and **M16** or **KSOM** medium for culture inside a CO<sub>2</sub> incubator. Many labs prefer the use of KSOM as it supports the development of embryos from many different mouse strains. All media can be prepared in-house according to the previously published method (Behringer, R., *et al.* Manipulating the mouse embryo: a laboratory manual, Fourth edition. ed. Cold Spring Harbor, New York: Cold Spring Harbor Laboratory Press; 2014.) or commercially purchased. Mouse oocytes and pre-implantation embryos are incubated in embryo culture medium pre-equilibrated for at least 4h in a CO<sub>2</sub> incubator (5% CO<sub>2</sub>, 37°C, 95% humidity).

1. 2 males of the same line should be sacrificed and sperm of better quality (based on a visual assessment after 60 min preincubation) should be utilized in an IVF procedure. If sperm of both males show the same quality both can be used in an IVF. Sacrifice two males and dissect the cauda epididymides. Transfer the 2 cauda epididymides (after removal of fat and blood) of each male to separate dishes into the oil next to the c-TYH (or FERTIUP® PM) drop. After nicking the tissue drag sperm with watchmaker forceps into the drop.
2. Sperm should be allowed capacitation for 60 min in c-TYH (or FERTIUP® PM);
3. 0.25 mM GSH concentration should be used in the mRVF. For this purpose:
  - Prepare a 100x  $\text{CaCl}_2 \times 2\text{H}_2\text{O}$  stock solution (310 mM) by dissolving 0.4558 g of  $\text{CaCl}_2 \times 2\text{H}_2\text{O}$  in 10 ml of embryo-grade water
  - Filter the solution through a 0.22  $\mu\text{m}$  filter and store aliquots at  $-20^\circ\text{C}$  for a maximum of 6 months
  - On the day of IVF thaw an aliquot of  $\text{CaCl}_2$  at room temperature
  - Add 150  $\mu\text{l}$  of 100x  $\text{CaCl}_2$  to 15 ml of RVF medium and mix gently
  - Place 1 ml of RVF medium supplemented with  $\text{CaCl}_2$  in a tube containing 30.7 mg of GSH and vortex
  - Add 10  $\mu\text{l}$  of this solution to 4 ml of RVF medium supplemented with  $\text{CaCl}_2$ , mix gently and filter using 0.22  $\mu\text{m}$  syringe end filter
4. Oocytes from a maximum of 5 females should be placed in a 200  $\mu\text{l}$  drop of mRVF;
5. 5  $\mu\text{l}$  of sperm suspension should be added to the oocytes. If the removal of cumulus cells assessed after 20 min of incubation is poor, indicating insufficient motility or concentration of sperm, additional 5  $\mu\text{l}$  of the sperm suspension should be transferred to the fertilization medium.

Aldrich Catalog #G8540

Aldrich Catalog #W1503

Aldrich Catalog #R7630

Dickinson Catalog #232100

**Bio Catalog #KYD-002-EX-X5**

Roth Catalog #671941

Aldrich Catalog #S5886

Aldrich Catalog #P5405

Aldrich Catalog #M2773

Aldrich Catalog #P5655

Aldrich Catalog #G6152

Aldrich Catalog #P7794

Aldrich Catalog #S1277

Aldrich Catalog #S8761

Scientific Catalog #11-360-070

Aldrich Catalog #C4555

Aldrich Catalog #P8136

Aldrich Catalog #C7902

Medical Catalog #K-RVFE-50

Aldrich Catalog #G4251

clear Minitube Catalog #13407/0010

Dish Corning Catalog #353001



- 8 Place the dish on a 37°C hot plate for 3 min and gently swirl every min for 20 sec to help the sperm disperse from the tissue
- 9 Divide the sperm suspension into 20 aliquots of 10 µl on a 10-cm culture dish lid avoiding carryover of paraffin oil into the aliquots (clean the pipette tip from the outside with a tissue to remove the oil each time before placing a 10 µl aliquot on the dish lid)
- 10 Aspirate each 10 µl drop into a separate freezing straw followed by 2.3 cm air
- 11 Seal the straws (e.g. with a heat sealer or metal balls) and place them in liquid nitrogen gas phase for 10 min. Note: We use a custom-made metal inlay for this purpose but self-made or purchased freezing canisters can be used, too (e.g. <http://card.medic.kumamoto-u.ac.jp/card/english/sigen/manual/spfreeze.html#canister> or KYD-S018 from Cosmo Bio)
- 12 Transfer the straws to the liquid nitrogen tank for long-term storage

#### Oocyte isolation

- 13 Ideally, the entire IVF procedure (oocyte preincubation, sperm preincubation and fertilization) should be performed in an incubator (5% CO<sub>2</sub>, 37°C) at 5% O<sub>2</sub> but atmospheric O<sub>2</sub> concentration have been shown to work well, too
- 14 Prepare a 35-cm culture dish with a 90 µl drop of mRVF (for preparation see Guidelines & Warnings) covered with oil (for oocytes from a maximum of 3 females) and equilibrate it for at least 20 min in an atmosphere of mixed gas (5% CO<sub>2</sub>, 5% O<sub>2</sub>, 37°C)
- 15 Collect oviducts from superovulated females 15 hours after the hCG injection and clean them in DPBS
- 16 Transfer the oviducts into the paraffin oil next to the 90 µl drop of mRVF
- 17 Release oocyte clutches into the oil by ripping the ampulla with forceps and drag them through the oil into the fertilization drop
- 18 Incubate oocytes for 50 min before adding the sperm suspension (at least 30 min and no longer than 60 min) in an atmosphere of mixed gas (5% CO<sub>2</sub>, 5% O<sub>2</sub>, 37°C)

#### Sperm thawing and capacitation

- 19 Prepare a 35-mm culture dish (for each IVF experiment) with 90 µl c-TYH drop (for preparation see Guidelines & Warnings) covered with paraffin oil and equilibrate it overnight in a CO<sub>2</sub> incubator (5% CO<sub>2</sub>, 37°C, 95% humidity) or at least 20 min in the morning of the IVF in an atmosphere of mixed gas (5% CO<sub>2</sub>, 5% O<sub>2</sub>, 37°C)

- 20 Remove the required straw(s) from long-term storage in liquid nitrogen on the day of IVF, place in a dewar with liquid nitrogen and then quickly transfer into a 37°C water bath for 10 min
- 21 Dry the straw(s) with a tissue and cut the sealed end and the labeled end of the straw below the cotton plug
- 22 Expel 10 µl sperm suspension into the center of a 90 µl c-TYH drop using a 1 ml syringe
- 23 Preincubate for 30 min in an atmosphere of mixed gas (5% CO<sub>2</sub>, 5% O<sub>2</sub>, 37°C) before the IVF procedure to allow capacitation of the sperm

#### *In vitro* fertilization

- 24 Add 10 µl of the sperm suspension taken from the edge of the c-TYH drop to the oocyte clutches with the help of a 200 µl cell-saver tip and incubate for 4 hours (at least 3 hours and no longer than 5 hours) in an atmosphere of mixed gas (5% CO<sub>2</sub>, 5% O<sub>2</sub>, 37°C)
- 25 Add another 10 µl of the sperm suspension to the fertilization medium if the removal of cumulus cells assessed after 20 min of incubation is poor, indicating insufficient motility or concentration of sperm

#### Embryo culture and transfer

- 26 Wash embryos after the IVF procedure through 10 drops of preincubated embryo culture medium (e.g. M16 or KSOM) and incubate overnight in embryo culture medium in groups of 15-50 embryos per drop (a 30 µl drop of embryo culture medium covered with paraffin oil) in a CO<sub>2</sub> incubator (5% CO<sub>2</sub>, 37°C, 95% humidity)
- 27 The day after insemination determine fertilization rates (a percentage of the total number of inseminated oocytes that developed to the 2-cell stage)
- 28 Transfer 2-cell embryos into the oviducts of pseudo-pregnant 0.5 dpc females