

Aug 29, 2024

Purification of FIP200 CTR-GFP

DOI

dx.doi.org/10.17504/protocols.io.j8nlk8866l5r/v1

Elias Adriaenssens¹, Eleonora Turco¹

¹Sascha Martens lab, University of Vienna, Max Perutz Labs - Vienna



Elias Adriaenssens

Sascha Martens lab, University of Vienna, Max Perutz Labs - ...

OPEN  ACCESS



DOI: dx.doi.org/10.17504/protocols.io.j8nlk8866l5r/v1

Protocol Citation: Elias Adriaenssens, Eleonora Turco 2024. Purification of FIP200 CTR-GFP. **protocols.io**

<https://dx.doi.org/10.17504/protocols.io.j8nlk8866l5r/v1>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: May 23, 2024

Last Modified: August 29, 2024

Protocol Integer ID: 101125

Keywords: ASAPCRN

Funders Acknowledgement:

**Aligning Science Across
Parkinson's (ASAP)**

Grant ID: ASAP-000350


**Marie Skłodowska-Curie
MSCA Postdoctoral
fellowship**

Grant ID: 101062916

Abstract

This protocol details the purification GFP-FIP200 (CTR).

Materials

 Rosetta™(DE3)pLysS Competent Cells - Novagen **Merck Catalog #70956-4**

Lysis buffer:

A	B
Tris HCl, pH 7.4	50 mM
NaCl	300 mM
MgCl ₂	2 mM
Glycerol	5%
Imidazole	10 mM
β-mercaptoethanol	2 mM

Wash buffer:

A	B
Tris HCl, pH 7.4	50 mM
NaCl	300 mM
Imidazole	10 mM
β-mercaptoethanol	2 mM

SEC buffer:

A	B
HEPES, pH 7.4	25 mM
NaCl	150 mM
DTT	1 mM



Purification procedure

1d 1h 45m 30s

- 1 To purify GFP-FIP200(CTR), as described previously (Turco et al. 2019 Mol Cell, PMID: 30853400), fuse the C-terminal domain of FIP200 (1458-1594aa) to a N-terminal 6xHis-TEV-GFP-tag through cloning into a pET-DUET1 vector (available on Addgene).
- 2 After the transformation of the pET-DUET1 vector encoding 6xHis-TEV-GFP-FIP200(CTR) in *E. coli* Rosetta pLysS cells (Novagen Cat# 70956-4), grow cells in 2x Tryptone Yeast extract (TY) medium at 37 °C until an OD₆₀₀ of 0.4 and then continued at 18 °C .
- 3 Once the cells reached an OD₆₀₀ of 0.8, induce the protein expression with 100 micromolar (μM) isopropyl β-D-1-thiogalactopyranoside (IPTG) for 16:00:00 at 18 °C . 16h
- 4 Collect the cells by centrifugation and resuspend in lysis buffer, complete EDTA-free protease inhibitors (Roche), CIP protease inhibitor (Sigma), and DNase (Sigma)).

Lysis buffer:

A	B
Tris HCl, pH 7.4	50 mM
NaCl	300 mM
MgCl ₂	2 mM
Glycerol	5%
Imidazole	10 mM
β-mercaptoethanol	2 mM

- 5 Sonicate cell lysates twice for 00:00:30 . 30s
- 6 Clear lysates by centrifugation at 18000 rpm, 4°C, 00:45:00 in a SORVAL RC6+ centrifuge with an F21S-8x50Y rotor (Thermo Scientific). 45m
- 7 Filter the supernatant through an 0.45 μm filter and loaded onto a pre-equilibrated 5 ml His-Trap HP column (Cytiva).



- 8 After binding of His tagged proteins to the column, wash the column with three column volumes of wash buffer.

**Wash buffer:**

A	B
Tris HCl, pH 7.4	50 mM
NaCl	300 mM
Imidazole	10 mM
β-mercaptoethanol	2 mM

- 9 Elute the proteins with a stepwise imidazole gradient (30, 75, 100, 150, 225, 300 mM).

- 10 Pool and incubate the fractions containing the 6xHis-TEV-GFP-FIP200(CTR) with TEV protease at 4 °C .

🕒 Overnight

8h



- 11 After cleave off the 6xHis tag , recapture 6xHis tag and His-tagged TEV protease with nickel beads for 01:00:00 at 4 °C .

1h

- 12 Pellet the beads by centrifugation and the supernatant, concentrate containing the GFP-FIP200(CTR) protein using a 30 kDa cut-off Amicon filter (Merck Millipore) and load onto a pre-equilibrated Superdex 200 Increase 10/300 GL column (Cytiva).



- 13 Elute the proteins with SEC buffer.

SEC buffer:

A	B
HEPES, pH 7.4	25 mM
NaCl	150 mM
DTT	1 mM


- 14 Analyse the fractions by SDS-PAGE and Coomassie staining.



- 15 Pool fractions containing purified GFP-FIP200(CTR).



16 After concentrating the purified protein, aliquot the protein and snap-frozen in liquid nitrogen.

17 Store proteins at  -80 °C .