



Feb 15, 2020

Quantitative analysis method of mRNA expression using dual RNA *in situ* hybridization by comparison with housekeeping genes

PLOS One

Byung Joon Seung¹¹Konkuk University

1 Works for me

dx.doi.org/10.17504/protocols.io.badcia2w Byung Joon Seung 

ABSTRACT

In retrospective formalin-fixed paraffin-embedded (FFPE) samples, RNA quality of FFPE tissues showed variation following the storage period and fixation process. So we suggest quantitative analysis by dual detection of target mRNA and housekeeping genes in one section and calculated the target gene/housekeeping gene ratio by two open-source image analysis programs. This assay potentially allows for reliable quantification of mRNA expression levels in retrospective FFPE samples.

EXTERNAL LINK

<https://doi.org/10.1371/journal.pone.0229031>

THIS PROTOCOL ACCOMPANIES THE FOLLOWING PUBLICATION

Seung B, Cho S, Kim S, Lim H, Sur J (2020) Quantitative analysis of *HER2* mRNA expression by RNA *in situ* hybridization in canine mammary gland tumors: Comparison with immunohistochemistry analysis. PLoS ONE 15(2): e0229031. doi: [10.1371/journal.pone.0229031](https://doi.org/10.1371/journal.pone.0229031)

GUIDELINES

RNA-ISH method was performed manually on FFPE tissues using the RNAscope[®] duplex assay (Advanced Cell Diagnostics, Hayward, CA, USA) according to the manufacture's instructions.

MATERIALS

NAME ▾


CATALOG # ▾

VENDOR ▾

RNAscope[®] duplex assay

RNA-ISH original image

- 1 RNA-ISH images were acquired at $\times 400$ magnification.




Example of zoomed image of RNA-ISH duplex image ($\times 400$ magnification)

Counting the blue dots

- 2 The original images were converted using Fiji program [1].


[1] Schindelin J, Arganda-Carreras I, Frise E, Kaynig V, Longair M, Pietzsch T, et al. Fiji: an open-source platform for biological-image analysis. Nat Methods. 2012;9(7):676-82.

- 3 Original images were converted using "Image→Color→Dichromacy→Tritanope" by Fiji program
Save as Tiff file.



- 4 Select regions of interest (ROI) in Tritanope filtered images.

To measure blue dots, dark spot detection mode in the spot detector of ICY program [2] was used.



Size of spots to

detect (scale 2).

Find adequate filter conditions for your research.


[2] De Chaumont F, Dallongeville S, Chenouard N, Hervé N, Pop S, Provoost T, et al. Icy: an open bioimage informatics platform for extended reproducible research. Nat Methods. 2012;9(7):690-6.

Counting the red dots

- 5 Tritanope-filtered images were converted to a CIELAB (RGB to CIELAB) and 'a' channel images were acquired as jpg files using Fiji program.

- 6 Select regions of interest (ROI) in 'a' channel images.

To measure red dots, white spot detection mode in the spot detector of ICY program was used.



Size of spots to detect (scale 3).

Find adequate filter conditions for your research.



This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited