

Jun 23, 2024

## Installing FIJI and SynBot (Windows Version)

 In 1 collection

DOI

**[dx.doi.org/10.17504/protocols.io.kqdg32qbpv25/v1](https://dx.doi.org/10.17504/protocols.io.kqdg32qbpv25/v1)**

Justin T Savage<sup>1</sup>

<sup>1</sup>Duke University

ASAP Collaborative Rese...

Eroglu\_Lab



Justin T Savage

Duke University

OPEN  ACCESS



DOI: **[dx.doi.org/10.17504/protocols.io.kqdg32qbpv25/v1](https://dx.doi.org/10.17504/protocols.io.kqdg32qbpv25/v1)**

**Protocol Citation:** Justin T Savage 2024. Installing FIJI and SynBot (Windows Version). **protocols.io**  
**<https://dx.doi.org/10.17504/protocols.io.kqdg32qbpv25/v1>**

**License:** This is an open access protocol distributed under the terms of the **[Creative Commons Attribution License](#)**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working

**Created:** June 12, 2024

**Last Modified:** June 23, 2024

**Protocol Integer ID:** 101692

## Abstract

Video instructions for installing FIJI and SynBot and a simple SynBot run for Windows operating system.



1

[https://www.youtube.com/embed/VZqD8DS\\_VfA?si=eJQ3iLecxYiwR4vR](https://www.youtube.com/embed/VZqD8DS_VfA?si=eJQ3iLecxYiwR4vR)

- 1.1 If video quality is poor, try watching directly from Youtube at [https://youtu.be/VZqD8DS\\_VfA?si=DFcoeyiLyKluYUUX](https://youtu.be/VZqD8DS_VfA?si=DFcoeyiLyKluYUUX)

## Install FIJI

- 2 0:00 - 0:40 Install FIJI from <https://fiji.sc> .

## Install SynBot

- 3 0:40 - 1:37 Download code from the SynBot GitHub repository [https://github.com/Eroglu-Lab/Syn\\_Bot](https://github.com/Eroglu-Lab/Syn_Bot) .
- 4 1:37 - 2:09 Move SynBot java plugins (JAR files) to the FIJI plugins folder. The **ilastik4ij\_Syn\_Bot-1.8.2-SNAPSHOT.jar** and **SynQuantExtra-1.2.9.jar** files are required for running SynBot.

## Run SynBot

- 5 2:09 - 2:22 Launch FIJI and open SynBot by dragging the SynBot.ijm file into FIJI.
- 6 2:22 - 2:49 Prepare your image folder. Create a folder for your experiment with subfolders for experimental groups. Each experimental group folder should contain only image files. Multiple experimental groups are not required, but your files must follow the Experiment/Group/Image structure.
- 7 2:49 - 3:41 Run SynBot (see the **SynBot manual thresholding** protocol for more information).
- 8 3:41 - 4:12 View results in Summary.csv file.