


Jun 17, 2021

Protocol for preparing minced *Artemia* nauplii

Ranjith Lakshmanan¹, Linga Prabu Dhanasekaran¹, Kalidas C¹, Mathanbabu A¹, Saravanan Raju¹¹ICAR-Central Marine Fisheries Research Institute

Ranjith Lakshmanan: Tuticorin Regional Station, ICAR-Central Marine Fisheries Research Institute

2 Works for me



Share

dx.doi.org/10.17504/protocols.io.bvirm4d6

ICAR TRS of CMFRI

Ranjith Lakshmanan

ICAR-Central Marine Fisheries Research Institute

ABSTRACT

Protocol for preparing minced *Artemia* nauplii as feed for young medusa of upside-down jellyfish

Cassiopea andromeda

- Laboratory hatched *Artemia* nauplii (Instar I & II stage) were filtered using a 250-micron sieve and concentrated to **15 mL** volume.
- The *Artemia* nauplii concentrate were centrifuging @ 8000 rpm for 3 minutes.
- The supernatant was removed and *Artemia* nauplii biomass was macerated using a sterile plastic pestle.
- Completely minced *Artemia* biomass was homogenously mixed with 5 mL of distilled water to make the concentration of 900-1100 nauplii/mL.
- Two drops (**100 µl**) of minced *Artemia* biomass was fed to young medusa.

DOI

dx.doi.org/10.17504/protocols.io.bvirm4d6

COLLECTION CITATION

Ranjith Lakshmanan, Linga Prabu Dhanasekaran, Kalidas C, Mathanbabu A, Saravanan Raju 2021.
Protocol for preparing minced *Artemia* nauplii. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.bvirm4d6>



KEYWORDS

Artemia, Jellyfish, Upside-down jellyfish, Feed, *Cassiopea andromeda*

LICENSE

This is an open access collection distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

Jun 04, 2021

LAST MODIFIED

Jun 17, 2021

COLLECTION INTEGER ID



50481

MATERIALS TEXT

1. *Artemia* nauplii (Instar I & II stage)
2. 250-micron mesh
3. Centrifuging machine (Lark)
4. Centrifuge tube (15 ml)
5. Glass pestle
6. Dropper

ABSTRACT

Protocol for preparing minced *Artemia* nauplii as feed for young medusa of upside-down jellyfish *Cassiopea andromeda*

- Laboratory hatched *Artemia* nauplii (Instar I & II stage) were filtered using a 250-micron sieve and concentrated to  **15 mL** volume.
- The *Artemia* nauplii concentrate were centrifuging @ 8000 rpm for 3 minutes.
- The supernatant was removed and *Artemia* nauplii biomass was macerated using a sterile plastic pestle.
- Completely minced *Artemia* biomass was homogenously mixed with 5 mL of distilled water to make the concentration of 900-1100 nauplii/mL.
- Two drops ( **100 µl**) of minced *Artemia* biomass was fed to young medusa.

FILES