



Jun 28, 2022

Ampicillin 100 mg/mL Stock Solution

Bailey Clark¹¹Florida State University iGEM

1 Works for me



Share

dx.doi.org/10.17504/protocols.io.ewov1no3ogr2/v1

FSU iGEM



Bec18

ABSTRACT

The abstract will be added later.

DOI

dx.doi.org/10.17504/protocols.io.ewov1no3ogr2/v1

PROTOCOL CITATION

Bailey Clark 2022. Ampicillin 100 mg/mL Stock Solution . **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.ewov1no3ogr2/v1>



LICENSE

————— This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

Feb 18, 2022

LAST MODIFIED

Jun 28, 2022


PROTOCOL INTEGER ID

58445

Stock Solution Preparation

5s

- 1 Weigh **100 mg** (**0.1 g**) of Ampicillin.

- 2 Add the Ampicillin to a 1.5 mL tube.
- 3 Add deionized water to the 1 mL level on the tube.
- 4 Vortex the tube for  00:00:05 , repeat 3 times.
- 5 Store the tube in the -20°C freezer.

5s