



May 17, 2022

# Mehlich 2 Solution

**Hsin-Mao Wu**<sup>1</sup><sup>1</sup>Kaohsiung Medical College[dx.doi.org/10.17504/protocols.io.q26g74xkkgwz/v1](https://dx.doi.org/10.17504/protocols.io.q26g74xkkgwz/v1) **Hsin-Mao Wu**

## DISCLAIMER – FOR INFORMATIONAL PURPOSES ONLY; USE AT YOUR OWN RISK

The protocol content here is for informational purposes only and does not constitute legal, medical, clinical, or safety advice, or otherwise; content added to [protocols.io](https://protocols.io) is not peer reviewed and may not have undergone a formal approval of any kind. Information presented in this protocol should not substitute for independent professional judgment, advice, diagnosis, or treatment. Any action you take or refrain from taking using or relying upon the information presented here is strictly at your own risk. You agree that neither the Company nor any of the authors, contributors, administrators, or anyone else associated with [protocols.io](https://protocols.io), can be held responsible for your use of the information contained in or linked to this protocol or any of our Sites/Apps and Services.

Soil extraction

DOI

[dx.doi.org/10.17504/protocols.io.q26g74xkkgwz/v1](https://dx.doi.org/10.17504/protocols.io.q26g74xkkgwz/v1)

Hsin-Mao Wu 2022. Mehlich 2 Solution. **protocols.io**  
<https://dx.doi.org/10.17504/protocols.io.q26g74xkkgwz/v1>



soil extraction

 protocol ,

Apr 27, 2022

May 17, 2022

61531

DISCLAIMER – FOR INFORMATIONAL PURPOSES ONLY; USE AT YOUR OWN RISK

The protocol content here is for informational purposes only and does not constitute legal, medical, clinical, or safety advice, or otherwise; content added to [protocols.io](https://protocols.io) is not peer reviewed and may not have undergone a formal approval of any kind. Information presented in this protocol should not substitute for independent professional judgment, advice, diagnosis, or treatment. Any action you take or refrain from taking using or relying upon the information presented here is strictly at your own risk. You agree that neither the Company nor any of the authors, contributors, administrators, or anyone else associated with [protocols.io](https://protocols.io), can be held responsible for your use of the information contained in or linked to this protocol or any of our Sites/Apps and Services.

Mehlich 2 Solution

1

在 **500 mL** DI Water中加入 **0.56 g**  $\text{NH}_4\text{F}$ 、**10.7 g**  $\text{NH}_4\text{Cl}$

2

再加 **1 mL**  $\text{HCl}$ 、**11.5 mL** 冰醋酸並加水至 **1 L**

此時pH值應為2.5左右